

Pivotal Therapy System™

The Pivotal Therapy System utilizes the progressively increased use of time and gravity to help reposition the spine and its soft tissue components. The components of this system help improve functional mobility and strengthen muscles which can reduce pain and improve spinal posture. Individuals should see improvement in pain, mobility and function within 3 to 4 weeks of use.

Guidelines For Use

Simply lie on the pivots in a pain free position. Specific guidelines regarding daily use of the Pivotal Therapy System should be provided by your healthcare professional.



Fabrication Enterprises®
250 Clearbrook Rd, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
FabEnt.com



MedEnvoy
NL-IM-000000248
Prinses Margrietplantsoen 33, Suite 123
2595 AM, The Hague (NL)



AJW Technology Consulting GmbH
Breite Strasse 3
40213 Düsseldorf (Germany)



Magister M™



Made in United States

© 2022 FEI, all rights reserved.

Pivotal Therapy System, Occiput, Cervical Pivot, Short Thoracic Pivot, Standard Thoracic Pivot, Lumbar Pivot, Magister, FEI and Fabrication Enterprises are trademarks of Goldberg. rev 0722



Occipivot™ (00-4502)

The Occipivot is a useful tool to reduce tension and stress at the base of the skull. It's beneficial for individuals with postural syndromes, myofascial pain, TMJ dysfunction, and headaches.

Instructions:

- Place the flat side of Occipivot on the floor, treatment table or bed.
- Lie on your back and position your neck against the concaved part of the Occipivot. Then rest the base of your skull on top of the fingertip nubs so they are in contact with your occipital muscles.
- Upon your first use, lie on the Occipivot in a pain-free position for no more than a couple minutes to determine tolerance. You can chose to move your head in flexion/extension, side bending, or rotation in attempt to find the pain-free position.
- Gradually increase usage time as recommended by your healthcare professional.

