



FingerWeightsTM



Our Mission

To help people improve finger and hand function across the world - We are FingerWeights™



FAQ's

Q. How do FingerWeights work?

A. FingerWeights™ utilize the accepted principles of progressive weight training practiced by athletes and prescribed by doctors to improve performance, prevent injury and rehabilitate.

Q. How long until I see results?

A. Feel the effects of FingerWeights after first use. Individual results will vary as with any exercise program.

Q. How often should I use the FingerWeights?

A. FingerWeights should be used every other day or every third day to give the muscles time to rebuild.

Q. Will FingerWeights make my fingers bigger?

A. No. The extensor and flexor muscles that open and close the fingers are located in the hands.

Q. How are FingerWeights different than other hand weights?

A. FingerWeights may be used while performing the targeted activity such as practicing guitar or playing video games. Other hand weight systems require a weighting object to be held making them impossible to use during practice or play.



Instructions

1

5



INSTRUCTIONS


www.fingerweights.com



Weight Rods

PARTS





FINGER MAP



A & B RING POSITIONS



WEIGHT CONFIGURATIONS

2

6







SIX-LEVEL TRAINING CYCLE

Designed for use with any practice or training regimen

- For Levels 1, 2, and 3 strap ring in Position A. Incease one level every 2 4 weeks.
- For Levels 4, 5, and 6 strap ring in Position B. Increase one level every 2 4 weeks.

LEVEL	DEVICE POSITION	WEIGHT CONFIG.	LEVEL DESCRIPTION
Level 1	6 6 A A	එ	Position A / 10 grams
Level 2	8890	Ô	Position A / 20 grams
Level 3	1 /	Ô	Position A / 30 grams
Level 4	B	එ	Position B / 10 grams
Level 5	Min	Ô	Position B / 20 grams
Level 6	' '	Ô	Position B / 30 grams

- After completing Level 6, decrease one level every 2 4 weeks.

ASSEMBLY









- Open strap guard.
- Hook *main eyelet* over either knob extending from ring.



Place *ring* on finger in Position A with capped end of ring facing fingertip.



 Wrap strap securely around finger to desired tension by hooking one of the *sizing* eyelets over opposite knob.



. Note: The knobs are designed to allow straps to overlap. Secur





IMPORTANT: For added safety and aerodynamics, capped end of weight chambers should always face fingertips.

REMOVING WEIGHTS

- With ring in hand (not attached to finger), insert weight ejector into small hole located in capped end of ring/weight chamber.
- · Push firmly until weight releases from chamber.





7

Close strap guard until it "clicks" into locked position.

1

Exercises

9

13

FINGER MAP



DEVICE POSITIONING



8 WEIGHT CONFIGURATIONS





WARNING: Choking Hazard -- Small parts. Not for children under 3 years.

SIX-LEVEL TRAINING SYSTEM

Designed for use with any practice or training regimen

- For Levels 1, 2, and 3 place ring in Position A. Increase one level every 4 weeks.
 For Levels 4, 5, and 6 place ring in Position B. Increase one level every 4 weeks.

LEVEL	DEVICE POSITION	WEIGHT CONFIG.	LEVEL DESCRIPTION
Level 1	0 , A A	Ô	Position A / 1 weight
Level 2	100	Ô	Position A / 2 weights
Level 3) /	Ô	Position A / 3 weights
Level 4	®	පී	Position B / 1 weight
Level 5	Min	Ô	Position B / 2 weights
Level 6	' /	Ô	Position B / 3 weights


- After completing Level 6, decrease one level every 4 weeks.
- Repeat entire cycle

10 HAND EXERCISE PROGRAM

Use Six-Level Training System in conjunction with this program

1. Wrist Curl

- Start with fingers relaxed
- Curl hands upward to full range of motion
- Hold for 3 seconds



- Reverse and curl downward
- Hold for 3 seconds
- Repeat 10 15 times

2. Fist To Claw

- Start with clenched fists
- Hold for 3 seconds





- Return to fist
- Repeat 10 15 times

12 3. Roof Top









- Return fingers to 90°
- Repeat 10 15 times

4. Finger Lift

- Start with hands on a flat surface fingers spaced evenly
- Lift each finger 5 10 times



· Reverse direction and repeat

5. Finger Walk

- Start with hands on a flat surface fingers together
- Beginning with the thumb, walk each finger inward
- Reverse direction and walk outward



6. Finger Circle

- Start with hands on a flat surface fingers spaced evenly
- Beginning with little fingers, lift and circle inward 5 10 times





14



Product Line

5-Finger Universal Set

Includes:

- 5 Finger Rings
- 15 Weights (3 per ring)
- Storage Pouch
- Weight Removal Tool
- Instructions for Assembly
- Exercise Program



Includes:

- 10 Finger Rings,
- 30 Weights (3 Per Ring)
- Storage Pouch
- Weight Removal Tool
- Instructions for Assembly
- Exercise Program





Therapy Pro



Tungsten Medical

