



Strength and Rehabilitation Products

Pedal Exerciser REF 10-0710

Adjustable resistance, light-weight pedal exerciser is pre-assembled.



10-0710



Portable lower and upper body exerciser. Use when exercise bicycle is not available or is contra-indicated. Wide leg spread and rubber non-marring feet keep exerciser stable even when patient has little muscle control. Wrap around adjustable footstraps for use with or without shoes.



Adjust the resistance level

Turn the tension knob to adjust the pedal resistance. Turn clockwise toward "High" to increase the resistance (makes pedaling more difficult) or counterclockwise towards "Low" to reduce the resistance (makes pedaling easier).



CAUTION Consult with your healthcare professional before beginning a new exercise regime with this product.



Fabrication Enterprises Inc.
250 Clearbrook Rd
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
www.FabEnt.com



Authorized CE representative:
AJW Technology Consulting GmbH
Königsallee 106
40215 Düsseldorf (Germany)