FLUIDOTHERAPY[®]

Dry Heat Therapy Units

User Manual

Models:

FLU110D

FLU110DE

FLU115D

FLU115DE

TABLE OF CONTENTS

Fluidotherapy® - Dry Heat Therapy

Foreword	1	Operation	24
About Dry Heat Therapy	2-5	Patient Preparation	24
Precautionary Instructions	2-4	Starting Treatment	24
Description of Device Markings	2	Stopping Treatment.....	24
Indications & Contraindications	5	Preventive Maintenance	25-32
Nomenclature	6-8	Daily Maintenance.....	25
FLU110D & FLU110DE Unit		Weekly Maintenance	26-29
Familiarization.....	6	ADDITIONAL FLUIDOTHERAPY PREVENTIVE	
FLU115D & FLU115DE Unit		MAINTENANCE REQUIREMENTS	
Familiarization.....	7	Quarterly	30
Operating Controls.....	8	Bi-Annual	30
Specifications	9-14	Annual.....	30
FLU110D & FLU110DE	9	As Needed	30
FLU115D & FLU115DE	10	Fluidotherapy Maintenance Record	31
Electromagnetic Compatibility Tables	11-14	Cleaning	32
Setup	15-23	Service	32
Treatment Mode Parameters.....	15-17	Accessories	33
Time Controlled Parameters	18	Replacement Accessories	33
Preference Mode Default Parameters.....	18-23	Warranty	34

FOREWORD

Fluidotherapy® - Dry Heat Therapy

Thank you for purchasing a Fluidotherapy Dry Heat Therapy Unit.

This manual contains general safety, operating, maintenance, and care instructions for the owners and operators of the Fluidotherapy Dry Heat Therapy Units.

The specifications put forth in this manual were in effect at the time of the publication. However, owing to DJO, LLC's policy of continuous improvement, changes to these specifications may be made at any time without obligation on the part of DJO, LLC.

Read, understand, and follow the information contained in this manual.

Stay current with the latest clinical developments in the field of Dry Heat Therapy and observe all applicable precautionary measures for treatment.

Keep informed on appropriate indications and contraindications for the use of Dry Heat Therapy.

This equipment is to be used only under the prescription and supervision of a licensed practitioner.

ABOUT DRY HEAT THERAPY

Fluidotherapy® - Dry Heat Therapy

PRECAUTIONARY INSTRUCTIONS

The precautionary instructions found in this section and throughout this manual are indicated by specific symbols. Understand these symbols and their definitions before operating this equipment. The definitions of these symbols are as follows:

= CAUTION - Text with a "CAUTION" indicator will explain possible safety infractions that could have the potential to cause minor to moderate injury or damage to equipment.

= WARNING - Text with a "WARNING" indicator will explain possible safety infractions that will potentially cause serious injury and equipment damage.

= DANGER - Text with a "DANGER" indicator will explain possible safety infractions that are imminently hazardous situations that would result in death or serious injury.

= EXPLOSION HAZARD - Text with an "Explosion Hazard" indicator will explain possible safety infractions if this equipment is used in the presence of flammable anesthetics.

= NON-IONIZING ELECTROMAGNETIC RADIATION - Text with a Non-Ionizing Electromagnetic Radiation" indicator informs the user of possible hazards resulting from elevated, potentially dangerous, levels of non-ionizing radiation.

= Protective Earth (ground)

NOTE: Helpful information marked as "NOTE" may be found throughout this manual. Each "NOTE" is beneficial in the aid of a particular area or function being described for this product.

DESCRIPTION OF DEVICE MARKINGS

Degree of Protection Against Electrical Shock

Complies with UL 2601-1 CSA C22.2 No. 601.1

Medical Electrical Equipment

Refer to Instruction Manual/Booklet

ABOUT DRY HEAT THERAPY

Fluidotherapy® - Dry Heat Therapy

CAUTION

- Read, understand, and practice the precautionary and operating instructions found in this manual. Know the limitations and hazards associated with using any electrical device. Observe the precautionary and operational decals placed on the unit.
- DO NOT operate this unit in an environment where other devices are being used that intentionally radiate electromagnetic energy in an unshielded manner. Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this manual.
- Portable and Mobile RF Communications equipment can affect Medical Electrical Equipment.
- This unit generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. Harmful interference to other devices can be determined by turning this unit on and off. Try to correct the interference using one or more of the following: reorient or relocate the receiving device, increase the separation between the equipment, connect the equipment to an outlet on a different circuit from that which the other device(s) are connected and consult the DJO, LLC Service Department for help.
- Clean Inlet Filter(s) daily before unit startup.
- Turn unit OFF before positioning a patient or removing a patient from the unit.
- Place the patient in a comfortable position allowing for correct placement of the limb being treated.
- Secure all entry ports before turning the unit ON.
- Use only Cellex® Dry Heat Medium in the Fluidotherapy units.
- Refill unit daily to proper fill level with Chattanooga Cellex Dry Heat Medium.

CAUTION

- Use only fingers to operate button controls on the control panel(s). Use of sharp objects such as pencils or pens will result in damage to the unit.
- Check unit temperature before treating patient to ensure correct temperature.
- This unit should be transported and stored in temperatures between -40 °F and 158 °F (-40 °C and -70 °C) with Relative Humidity between 10% - 100% to prevent damage to the unit or its components.
- This unit should be operated in temperatures between 110 °F and 125 °F (43 °C and 52 °C) to prevent damage to the unit or its components.
- DO NOT operate the unit when connected to any unit other than Chattanooga devices.
- Change Cellex Dry Heat Medium every six (6) months.
- Do not use accessories other than those supplied with the unit, or recommended by DJO, LLC. The safety of other products has not been established, and their use could result in injury to the patient.

ABOUT DRY HEAT THERAPY

Fluidotherapy® - Dry Heat Therapy

WARNING

- Federal law restricts this device to sale by, or on the order of, a physician or licensed practitioner. This device should be used only under the continued supervision of a physician or licensed practitioner.
- Make certain the unit is electrically grounded by connecting only to a grounded electrical service receptacle conforming to the applicable national and local electrical codes.
- Before administering any treatment to a patient, you should become acquainted with the operating procedures for each mode of treatment available, as well as the indications, contraindications, warnings, and precautions. Consult other resources for additional information regarding the application of Dry Heat Therapy.
- The Fluidotherapy Dry Heat Therapy Unit should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Fluidotherapy Dry Heat Therapy Unit should be observed to verify normal operation in the configuration in which it will be used.
- For continued protection against fire hazard, replace fuses only with ones of the same type and rating.
- This device should be kept away from children.
- To prevent electrical shock, disconnect the unit from the power source before attempting any maintenance procedures.
- Use only processed dry heat medium in the unit such as Cellex to prevent excessive dusting.
- Use only accessories that are specially designed for this unit. Do not use accessories manufactured by other companies on this unit. DJO, LLC is not responsible for any consequence resulting from using products manufactured by other companies. The use of other accessories or cables may result in increased emissions or decreased immunity of this unit.
- Dispose of all products in accordance with local and national regulations and codes.
- Fluidotherapy equipment not in use should be protected against unqualified use.

DANGER

- The Fluidotherapy Dry Heat Therapy Unit is not suitable in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.
- Perform all Preventive Maintenance as described in this manual. Failure to perform the Preventive Maintenance could result in the Cellex medium entering the heat chamber of the unit and cause severe injury to patients as well as smoke damage to the facility and the Fluidotherapy unit.
- The solvents of adhesives and flammable solutions used for cleaning and disinfecting the unit should be allowed to evaporate before the unit is used.

ABOUT DRY HEAT THERAPY

Fluidotherapy® - Dry Heat Therapy

Indications & Contraindications

Indications

- Relief of local pain
- Treatment of local blood flow insufficiency
- Treatment in range of motion when combined with exercise
- Treatment for symptoms of non-rheumatoid arthritis

Contraindications

- This device should not be used for symptomatic pain relief unless etiology is established or unless a pain syndrome has been diagnosed.
- This device should not be used when cancerous lesions or open wounds are present in the treatment area.
- Other counterindications are patients suspected of carrying serious infectious disease or disease where it is advisable, for general medical purposes, to suppress heat or fevers.


WARNING

- Adequate precautions should be taken when treating individuals with suspected or diagnosed medical conditions or diseases such as heart problems, epilepsy, diabetes, etc.
- Prior to treatment, consult a medical professional familiar with the precautionary measures to be taken for patients that may experience allergic reactions to dust and pollen.

NOMENCLATURE

Fluidotherapy® - Dry Heat Therapy

FLU110D & 110DE Unit Familiarization

1. RESERVOIR LID & TREATMENT VIEWING WINDOW

Access to add/change medium, view patient treatment, and clinician access to patient treatment area

2. MEDIUM/TREATMENT RESERVOIR

Treatment access and medium reservoir

3. HAND/ARM ACCESS PORT

Treatment port for Hand/Arm Treatment (end of unit)

4. REPLACEABLE INLET FILTERS

Filters room air entering the Fluidotherapy unit (1 each side)

5. FOUR POINT BASE

Casters included to convert to a Mobile Base

6. BLOWER HOUSING

Houses Blower & Heater

7. CONTROL PANEL

See [page 8](#) for detailed description of controls


8. ELBOW/FOOT TREATMENT ACCESS PORT

Treatment port for Elbow and Foot

NOMENCLATURE

Fluidotherapy® - Dry Heat Therapy

FLU115D & FLU115DE Unit Familiarization

1. ELBOW/FOOT ACCESS PORT

Treatment port for Elbow and Foot Treatment, access to add/change medium, clinician access to patient treatment area

2. HAND/ARM ACCESS PORTS

Treatment ports for Hand/Arm Treatment (end of unit)

3. INLET FILTERS

Filters room air entering the Fluidotherapy unit (1 each side)

4. MEDIUM/TREATMENT TUB

Treatment access and medium reservoir

5. MOBILE BASE

6. LOCKING CASTERS

7. BLOWER HOUSING

Houses Blower & Heater


8. CONTROL PANEL

See [page 8](#) for detailed description of controls

NOMENCLATURE

Fluidotherapy® - Dry Heat Therapy

Operating Controls (All Units)

1. DISPLAY

Displays Treatment Time, Temperature, Airspeed, Pulse Time, Preheat Timer, and Clock when the respective Mode is indicated

2. PREHEAT TIMER ON

Indicator light for Preheat Timer

3. INCREASE/HOUR

Use to increase Mode parameters and use to set the Hour when setting Clock

4. DECREASE/MINUTE

Use to decrease Mode parameters and use to set the Minutes when setting Clock

5. STOP

Use to stop treatment and Preheat Timer

6. MOTOR FUSE

7. START

Use to start treatment

8. INDICATOR BAR

Indicates Mode as they are chosen

9. HEATER FUSE

10. PREHEAT TIMER BUTTON

Press once to turn Preheat Timer On. Press twice to start Preheat function.

11. MODE SELECT

Use to select desired Mode

12. PM

PM indicator for Clock

SPECIFICATIONS

Fluidotherapy® - Dry Heat Therapy

FLU110D & FLU110DE

MODES OF OPERATION

Continuous (Default)

Variable Adjustments for Time, Temp, and Air Speed

Pulse Mode.....OFF to 6 Sec ON/OFF

TREATMENT TIME.....1 to 99 minutes

OPERATING TEMPERATURE.....110 °F (43.3 °C) to 125 °F (51.6 °C)

AIR SPEED.....0% to 100% (5% increments)

PREHEAT TIMER.....115 °F (46.1 °C) with 50% Air Flow

MEDIUM CAPACITY.....Approximately 30 lb (13.6 kg)

INPUT POWER (FLU110D).....120V, 50/60 Hz, 12A

INPUT POWER (FLU110DE).....230-240V, 50/60 Hz, 8A

PHYSICAL DIMENSIONS

Cabinet Length.....34.0" (86.4 cm)

Cabinet Width.....11.5" (29.2 cm)

Height.....33.0" (83.8 cm)

Weight.....70 lbs (31.7 kg) Less Medium

Shipping Weight.....100 lbs (45.4 kg)

DEGREE OF PROTECTION AGAINST INGRESS OF WATER.....IPX0

ELECTRICAL CLASS.....CLASS I

DEGREE OF PROTECTION AGAINST ELECTRICAL SHOCK.....TYPE B ⚡

POWER CORD.....2 meter cable length, 14 AWG, shielded

WARNING

Make certain the unit is electrically grounded by connecting only to a grounded electrical service receptacle conforming to the applicable national and local electrical codes.

SPECIFICATIONS

Fluidotherapy® - Dry Heat Therapy

FLU115D & FLU115DE

MODES OF OPERATION

Continuous (Default)

Variable Adjustments for Time, Temp, and Air Speed

Pulse Mode.....OFF to 6 Sec ON/OFF

TREATMENT TIME.....1 to 99 minutes

OPERATING TEMPERATURE.....110 °F (43.3 °C) to 125 °F (51.6 °C)

AIR SPEED.....0% to 100% (5% increments)

PREHEAT TIMER.....115 °F (46.1 °C) with 50% Air Flow

MEDIUM CAPACITY.....Approximately 40 lbs (18.1 kg)

INPUT POWER (FLU115D).....120V, 50/60 Hz, 12A

INPUT POWER (FLU115DE).....230-240V, 50/60 Hz, 10A

PHYSICAL DIMENSIONS

Cabinet Length.....34.0" (86.4 cm)

Cabinet Width.....18.5" (47.0 cm)

Height.....33.0" (83.8 cm)

Weight.....60 lbs (27.2kg) Less Medium

Shipping Weight.....100 lbs (45.4 kg)

DEGREE OF PROTECTION AGAINST INGRESS OF WATER.....IPX0

ELECTRICAL CLASS.....CLASS I

DEGREE OF PROTECTION AGAINST ELECTRICAL SHOCK.....TYPE B ⚡

POWER CORD.....2 meter cable length, 14 AWG, shielded

WARNING

Make certain the unit is electrically grounded by connecting only to a grounded electrical service receptacle conforming to the applicable national and local electrical codes.

SPECIFICATIONS

Fluidotherapy® - Dry Heat Therapy

ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

TABLE 1: GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

The Fluidotherapy - Dry Heat Therapy Unit is intended for use in the electromagnetic environment specified in the table below. The user of the Fluidotherapy - Dry Heat Therapy Unit should assure that it is used in such an environment.

Emission Tests	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The Fluidotherapy - Dry Heat Therapy Unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The Fluidotherapy - Dry Heat Therapy Unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations IEC 61000-3-3	Complies	

SPECIFICATIONS

Fluidotherapy® - Dry Heat Therapy

ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES (CONTINUED)

TABLE 2 & 3: GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The Fluidotherapy - Dry Heat Therapy Unit is intended for use in the electromagnetic environment specified in the table below. The user of the Fluidotherapy - Dry Heat Therapy Unit should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Fluidotherapy - Dry Heat Therapy Unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance:
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	$d = \frac{[3.5]\sqrt{P}}{V_1}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \frac{[3.5]\sqrt{P}}{E_1}$ 80 MHz to 800 MHz $d = \frac{[Z]\sqrt{P}}{E_1}$ 800 MHz to 2.5 GHz
			where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Fluidotherapy - Dry Heat Therapy Unit is used exceeds the applicable RF compliance level above, the Fluidotherapy - Dry Heat Therapy Unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Fluidotherapy - Dry Heat Therapy Unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

SPECIFICATIONS

Fluidotherapy® - Dry Heat Therapy

ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES (CONTINUED)

TABLE 2 & 3: GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The Fluidotherapy - Dry Heat Therapy Unit is intended for use in the electromagnetic environment specified in the table below. The user of the Fluidotherapy - Dry Heat Therapy Unit should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge ⁽¹⁾ IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Fluidotherapy - Dry Heat Therapy Unit requires continued operation during power mains interruptions, it is recommended that the Fluidotherapy - Dry Heat Therapy Unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

SPECIFICATIONS

Fluidotherapy® - Dry Heat Therapy

ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES (CONTINUED)

TABLE 4: RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE FLUIDOTHERAPY - DRY HEAT THERAPY UNIT (CONTINUED)

The Fluidotherapy - Dry Heat Therapy Unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Fluidotherapy - Dry Heat Therapy Unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Fluidotherapy - Dry Heat Therapy Unit as recommended below.

Rated Maximum Output Power of Transmitter W	Separation Distance According to Frequency of Transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \frac{[3.5]\sqrt{P}}{V_1}$	$d = \frac{[3.5]\sqrt{P}}{E_1}$	$d = \frac{[7]\sqrt{P}}{E_1}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

SETUP

Fluidotherapy® - Dry Heat Therapy

TREATMENT MODE PARAMETERS

ATTENTION!!!

**ATTRITION OF PARTICLES DURING SHIPMENT
MAY RESULT IN INITIAL DUSTING AND SMALL
PARTICLE RELEASE. DUST AND ODOR WILL
DISSIPATE WITH USE.**

Use the following instructions to set the various Mode Parameters to the desired settings.

TREATMENT TIME

Press the "PRESS TO CHANGE MODE" button until the indicator light is beside "TREATMENT TIME".

NOTE:

The default time of "20:00" will be displayed.

Press the "INCREASE" or "DECREASE" buttons until desired Treatment Time is displayed.

NOTE:

The "TREATMENT TIME" can be adjusted in one minute increments from 1 to 99.

TREATMENT TEMPERATURE

Press the "PRESS TO CHANGE MODE" button until the indicator light is beside "TEMPERATURE".

NOTE:

The default temperature of "115 °F (46.1 °C)" will be displayed when unit is first plugged into wall outlet.

The "F" or "C" and "TEMPERATURE" indicator will flash while the programmed Treatment Temperature is being displayed. When the Bed Temperature is displayed, the "F" or "C" and the "TEMPERATURE" indicator will illuminate steadily.

SETUP

Fluidotherapy® - Dry Heat Therapy

TREATMENT MODE PARAMETERS (CONTINUED)

TREATMENT TEMPERATURE (CONTINUED)

Press the "INCREASE" or "DECREASE" buttons until desired Treatment Temperature is displayed.

NOTE:

°F can be changed to °C and vice versa as desired. See [page 20](#) for instructions. Treatment Temperature can be adjusted in 1° increments from 88 °F to 130 °F (31.1 °C to 54 °C).

Press the "INCREASE" or "DECREASE" buttons until desired Air Speed is displayed.

NOTE:

The "AIR SPEED" is adjusted in increments of 5 from 0 to 100.

AIR SPEED

Press the "PRESS TO CHANGE MODE" button until the indicator light is beside "AIR SPEED".

NOTE:

The default speed of "50" will be displayed.

PULSE TIME

Press the "PRESS TO CHANGE MODE" button until the indicator light is beside "PULSE TIME".

NOTE:

The "PULSE TIME" allows the unit to operate by pulsing the medium during treatment. If selected and set, the medium will pulse on and off in equal increments. The Factory Default is "OFF".

EXAMPLE: "PULSE TIME" is set at "5". The unit will pulse the medium during treatment, five seconds on and five seconds off.

SETUP

Fluidotherapy® - Dry Heat Therapy

TREATMENT MODE PARAMETERS (CONTINUED)

PULSE TIME (CONTINUED)

Press the "INCREASE" or "DECREASE" buttons until desired Pulse Time is displayed.

NOTE:

"PULSE TIME" is adjusted in one second increments from "OFF" to 6 seconds.

STOP TREATMENT

Press the "STOP" button.

START TREATMENT

Press the "START" button.

SETUP

Fluidotherapy® - Dry Heat Therapy

TIME CONTROLLED PARAMETERS

PREHEAT TIMER

Press the "PRESS TO CHANGE MODE" button until the indicator light is beside "PREHEAT TIMER".

NOTE:

The "PREHEAT TIMER" allows the unit to preheat the medium. The unit will run at 50% Air Speed until the unit reaches the default preheat temperature or 90 minutes, whichever comes first. The "PREHEAT TIMER" will operate Monday through Friday only and can be set to automatically start at a predetermined time and heat the medium to a predetermined temperature. See [page 21](#) for setting the default parameters of the "PREHEAT TIMER".

Press the "PRESS ONCE TO TURN ON PREHEAT TIMER" button to illuminate the "PREHEAT TIMER ON" indicator light. Press again to start the "PREHEAT TIMER".

NOTE:

The "PREHEAT TIMER" indicator light must be set at the end of each day in order for it to automatically come on the next day. The "STOP" button is used to turn off the Preheat if desired.

CLOCK

Press the "PRESS TO CHANGE MODE" button until the indicator light is beside "CLOCK".

NOTE:

The "PREHEAT TIMER" will not function until the Clock is set.

Press the "INCREASE" button to set hours. Press the "DECREASE" button to set minutes.

NOTE:

The "PM" indicator light will illuminate when PM hours are reached.

SETUP

Fluidotherapy® - Dry Heat Therapy


PREFERENCE MODE DEFAULT PARAMETERS

ENTERING "PrEF" MODE (PREFERENCE MODE)

Simultaneously press the "PRESS TO CHANGE MODE", "INCREASE", and "DECREASE" buttons. "PrEF" will be displayed.

NOTE:

The "PrEF" mode allows the user to change the default settings.

TREATMENT TIME DEFAULT

Press the "PRESS TO CHANGE MODE" button until the indicator light is flashing beside "TREATMENT TIME".

Press the "INCREASE" and "DECREASE" buttons to set unit to the desired default treatment time.

Press the "PRESS TO CHANGE MODE" button to save the new setting.

NOTE:


The default setting is 20:00 minutes.

SETUP

Fluidotherapy® - Dry Heat Therapy

PREFERENCE MODE DEFAULT PARAMETERS (CONTINUED)

TEMPERATURE READING DEFAULT

Press the "PRESS TO CHANGE MODE" button until the indicator light is flashing beside "TEMPERATURE".

TREATMENT TEMPERATURE DEFAULT

Press the "PRESS TO CHANGE MODE" button until the indicator light is flashing beside "TEMPERATURE" and the existing default temperature is displayed.

Press the "INCREASE" button to set unit to display °F or °C as desired for the default setting.

Press the "PRESS TO CHANGE MODE" button to save the new setting.

NOTE:

The default setting is °F.

Press the "INCREASE" and "DECREASE" buttons to set unit to the desired default treatment temperature.

Operating Temperature Range is:
110 °F - 125 °F (43.3 °C - 51.6 °C)

Press the "PRESS TO CHANGE MODE" button to save the new setting.

SETUP

Fluidotherapy® - Dry Heat Therapy

PREFERENCE MODE DEFAULT PARAMETERS (CONTINUED)

PULSE MODE (Enable/Disable)

Press the "PRESS TO CHANGE MODE" button until the indicator light is flashing beside "PULSE TIME".

"EnAb" (enable) or "dISA" (disable) will display.

PREHEAT TIMER MODE (Enable/Disable)

Press the "PRESS TO CHANGE MODE" button until the indicator light is flashing beside "PREHEAT TIMER".

"EnAb" (enable) or "dISA" (disable) will display.

Press the "INCREASE" button to set unit to display "EnAb" or "dISA" as desired.

Press the "PRESS TO CHANGE MODE" button to save the new setting.

Press the "INCREASE" button to set unit to display "EnAb" or "dISA" as desired.

Press the "PRESS TO CHANGE MODE" button to save the new setting.

NOTE:

If Clock is disabled, the "PREHEAT TIMER" will automatically be disabled.

SETUP

Fluidotherapy® - Dry Heat Therapy

PREFERENCE MODE DEFAULT PARAMETERS (CONTINUED)

CLOCK (Enable/Disable)

Press the "PRESS TO CHANGE MODE" button until the indicator light is flashing beside "CLOCK".

"EnAb" (enable) or "dISA" (disable) will display.

Press the "INCREASE" button to set unit to display "EnAb" or "dISA" as desired.

Press the "PRESS TO CHANGE MODE" button to save the new setting.

NOTE:

If Clock is disabled, the "PREHEAT TIMER" will automatically be disabled.

SETTING DAY OF WEEK

Press the "PRESS TO CHANGE MODE" button until "dAy_" (_ = day code) is displayed.

DAY CODE CHART

WEEKDAY	DAY CODE
Sunday	dAyu
Monday	dAyn
Tuesday	dAyE
Wednesday	dAyD
Thursday	dAyr
Friday	dAyF
Saturday	dAyA

Press the "INCREASE" or "DECREASE" button to set unit to display the desired day code.

Press the "PRESS TO CHANGE MODE" button to save the new setting.

SETUP


Fluidotherapy® - Dry Heat Therapy

PREFERENCE MODE DEFAULT PARAMETERS (CONTINUED)

EXITING "PrEF" MODE

Press the "PRESS TO CHANGE MODE" buttons until "PrEF" is displayed.

Simultaneously press the "INCREASE", "DECREASE", and "STOP" buttons. The default treatment time will display and the indicator will be beside "TREATMENT TIME".

OPERATION

Fluidotherapy® - Dry Heat Therapy

PATIENT PREPARATION

- Clean the patient treatment area thoroughly with an antimicrobial soap and clean water as per industry, facility, regulatory standards, and universal skin washing procedures.
- Following the skin washing procedure, apply a hospital grade antiseptic skin cleanser, according to the cleanser manufacturer's recommended instructions for use.
- Reference the maintenance schedule and make certain the preventive maintenance has been performed prior to starting the unit.

STARTING TREATMENT

STARTING

Make certain all control panel settings have been made, unit is preheated, and patient is in proper position with sleeve(s) secure prior to starting the treatment. Refer to [pages 15 - 17](#) for setup of the unit. Press "START" button to begin treatment.

STOPPING TREATMENT

STOPPING

Treatment will automatically stop when the treatment time reaches zero. Should it be desired to stop treatment before the timer reaches zero, push the "STOP" button.

NOTE:

Should the treatment be stopped before the timer reaches zero, it will be necessary to re-set the "TREATMENT TIME" in order to complete the prescribed treatment.

PREVENTIVE MAINTENANCE

Fluidotherapy® - Dry Heat Therapy

DAILY MAINTENANCE


WARNING

Before any maintenance is performed or attempted, unplug the unit from the power source to prevent the possibility of electrical shock.

CLEAN INLET FILTERS

At the end of each work day, unplug the unit and clean the Inlet Filters on the unit. Carefully remove the filter retainer and wash the filter and screen with a mild antibacterial soap and water. Thoroughly dry the filter and screen before placing back on the unit.

NOTE:

Should your unit have the earlier style inlet filters, gently clean the filter using a soft bristle brush. Be careful not to puncture or damage the filter. Should the filter become damaged, torn, or punctured, call your dealer for replacement of the filter before resuming operation.

REFILL WITH CELLEX MEDIUM

Refill the unit with Cellex Dry Heat Medium to approximately one (1) inch above the bottom of the arm treatment ports of the unit.

NOTE:

Use of other than Cellex Medium may cause premature failure of the Fluidotherapy unit(s).

INSPECT SLEEVE CONDITION

Inspect the port sleeves for tears, rips, and weak seams. Replace all sleeves that show signs of tears, rips, weak, or loose seams, or excessive wear. Keeping the sleeves in excellent condition prevents excessive spillage of the Cellex medium and prevents the medium from entering the heat chamber of the unit.

ARM PREHEAT TIMER

Plug the unit back into an approved power source. If the "Preheat Timer" is used, press the "PREHEAT TIMER" button once and the "PREHEAT TIMER ON" indicator will illuminate.

NOTE:

The "PREHEAT TIMER ON" indicator must be illuminated in order for the Preheat Timer to automatically come on at the default time the next morning.


The Fluidotherapy Preheat Timer does not operate on Saturday or Sunday. However, for Monday morning preheat, set the unit as described above on Friday evening.

PREVENTIVE MAINTENANCE

Fluidotherapy® - Dry Heat Therapy

WEEKLY MAINTENANCE

Each week all sleeves of the Fluidotherapy unit should be laundered in a mild antibacterial detergent. Allow the sleeves to air dry or dry on a low temperature setting. Drying the sleeves in high temperatures could cause the sleeves to shrink or become distorted, resulting in the sleeve(s) not properly fitting when placing them back onto the unit.

REMOVING TOP SLEEVES (All Models)


Unplug the unit from the power source.

Carefully remove the rectangular bezel(s) from the top of the unit.

Carefully pull the hook and loop fasteners apart to remove the sleeve from the unit.

Replace in reverse order making certain the corners are well seated to form the seal required to prevent excessive escape of medium.

REMOVING END SLEEVE (Model FLU110D & FLU110DE)

Carefully pull the hook and loop fastener on the sleeve to remove sleeve from unit.

Replace in reverse order making certain the corners are well seated to form the seal required to prevent excessive escape of medium.

PREVENTIVE MAINTENANCE

Fluidotherapy® - Dry Heat Therapy

WEEKLY MAINTENANCE (CONTINUED)

REMOVING END SLEEVES (Model FLU115D & FLU115DE)

NOTE:


The sleeve in these photos has been modified (cut away) for clarity.

Remove the top sleeves to gain access to the inside of the unit. Refer to [page 26](#).

Pull the medium back and away from the arm ports.


NOTE:

It is extremely important to move the medium away from the ports to prevent the medium entering the space between the inner and outer walls of the unit when the sleeves are removed and replaced.

Dust away all Cellex medium from inside the sleeve.

Pull the retainer ring from the sleeve and lay aside for re-installation.

Pull the sleeve away from the sleeve extrusion in the inner housing of the unit.

NOTE:

There is a ring folded into the sleeve that fits over the sleeve extrusion on the inside of the unit.

From the outside of the unit, reach into the sleeve and remove the other sleeve retainer ring.

Cellex medium entering the heat chamber of the unit can cause severe injury to patients as well as smoke damage to the facility and the Fluidotherapy unit.

Carefully pull the sleeve and ring out of the unit, avoiding Cellex medium from entering the cavity between the two housings of the unit.

PREVENTIVE MAINTENANCE

Fluidotherapy® - Dry Heat Therapy

WEEKLY MAINTENANCE (CONTINUED)

Remove the ring from the folded end of the sleeve.
Launder sleeves as described on [page 26](#).

Vacuum out any of the Cellex medium that may have entered the cavity between the inner and outer walls of the unit.

WARNING


Removal of all Cellex medium that has entered the cavity is essential to maintain proper and safe operation of the unit. If the medium reaches the heating chamber of the unit, immediately call a certified technician for removal of any suspected medium in the heat chamber. Do not place the unit into service if medium is suspected to be in the heat chamber.

INSTALLING 115 Model SLEEVE

Fold the elastic end of the sleeve over the ring.


Place sleeve through port hole.

Place the ring/sleeve assembly over the sleeve extrusion on the inside of the unit.

NOTE:

Position the sleeve seam toward the top of the unit and push the ring/sleeve assembly completely onto the sleeve extrusion.

Install the inner sleeve retainer ring making certain the ring split will be toward the top of the unit.

Push the retainer ring over the ring/sleeve assembly until it is completely seated.


NOTE:

The retainer ring should fit tightly over the sleeve. If necessary, use a pair of pliers to tighten the ring in the area where it slides over the sleeve.

PREVENTIVE MAINTENANCE

Fluidotherapy® - Dry Heat Therapy

WEEKLY MAINTENANCE (CONTINUED)

From the outside of the unit, reach through the sleeve and install the outer sleeve retaining ring.

NOTE:

Start installation of the outer ring at the top of the unit. As the ring is being seated, pull the area of the sleeve tight. This will aid in the prevention of excessive medium spillage and patient comfort.

PREVENTIVE MAINTENANCE

Fluidotherapy® - Dry Heat Therapy

ADDITIONAL FLUIDOTHERAPY PREVENTIVE MAINTENANCE REQUIREMENTS

The following additional preventive maintenance requirements must be scheduled and performed as described to ensure that the unit is operating efficiently, safely, and functioning at optimum level. A blank Maintenance Record is provided on [page 31](#) to aid in the scheduling and record keeping of this prescribed preventive maintenance program. The following preventive maintenance procedures must be performed by a DJO, LLC qualified service technician trained in the maintenance requirements of the Chattanooga Fluidotherapy units.

QUARTERLY (Every 3 Months)

The following Preventive Maintenance must be performed on all Fluidotherapy units quarterly by a certified service technician.

- **INTERNAL CAVITY INSPECTION AND CLEANING**
- **FULL FUNCTIONAL AND PERFORMANCE TESTS**

BI-ANNUAL (Every 6 Months)

The following Preventive Maintenance must be performed on all Fluidotherapy units every six months in addition to the Quarterly Maintenance requirements by a certified service technician.

- **CHANGE CELLEX MEDIUM**

ANNUAL (Once per Year)

The following Preventive Maintenance must be performed on all Fluidotherapy units annually in addition to the Quarterly and Bi-Annual Maintenance requirements by a certified service technician.

- **CALIBRATION**

AS NEEDED

The following Preventive Maintenance must be performed only if performance test results indicate replacement is necessary.

- **INTAKE FILTER REPLACEMENT**
- **DISTRIBUTOR REPLACEMENT**
- **BLOWER MOTOR(S) REPLACEMENT**

DANGER

- If the level of Cellex medium suddenly drops one inch or more below the operating level of the unit, immediately pull the unit out of service and contact a certified service technician.
- A sudden level drop in the Cellex medium may indicate that medium is leaking into the heating chamber of the unit and must be repaired before the unit is placed back into service.

PREVENTIVE MAINTENANCE

Fluidotherapy® - Dry Heat Therapy

FLUIDOTHERAPY MAINTENANCE RECORD

UNIT SERIAL NUMBER _____ UNIT MODEL NUMBER _____ DATE PLACED IN SERVICE _____

DEALER _____ PHONE _____ CONTACT _____

DATE	MAINTENANCE PERFORMED	TECH INITIALS

DATE	MAINTENANCE PERFORMED	TECH INITIALS

Completing this form:

"DATE"-Date Service is performed

"MAINTENANCE PERFORMED"- Quarterly, Bi-Annual or Annual

"TECH INITIALS"- Certified Tech's Initials

PREVENTIVE MAINTENANCE

Fluidotherapy® - Dry Heat Therapy

CLEANING THE FLUIDOTHERAPY UNIT

After each use, clean the unit and its accessories using a soft, clean cloth dampened with water and a mild antibacterial detergent. Avoid the use of abrasive materials and cleaning solvents.

SERVICE

When the Fluidotherapy Dry Heat Therapy unit requires service or preventive maintenance, contact the selling dealer or DJO, LLC Service Department.

All units returned to the factory for service must include the following:

WARRANTY REPAIR/OUT OF WARRANTY REPAIR

1. Written statement containing the following information:
 - RA Number- Obtain from Factory
 - Unit Model Number
 - Unit Serial Number
 - Contact person with Phone and Fax Numbers
 - Billing Address (for Out of Warranty Repair)
 - Shipping Address (Where to Ship Unit after Repair)
 - Detailed Description of Problem or Symptoms
2. Copy of original invoice issued at purchase of the unit.
3. Ship unit to Factory in the original container with all accessories and information as required in item 1 above to:
DJO, LLC
Chattanooga Repair Center
47492 SD Hwy 22
PO Box 709
Clear Lake, SD 57226 USA
chattgroup.com

Service to these units should be performed only by service technicians certified by DJO, LLC.

Council Directive 2002/96/EC concerning Waste Electrical and Electronic Equipment (WEEE). Indicates a requirement not to dispose of WEEE as municipal waste. Contact your local distributor for information regarding disposal of the unit and accessories.

ACCESSORIES

Fluidotherapy® - Dry Heat Therapy

REPLACEMENT ACCESSORIES

When ordering additional accessories for Chattanooga Fluidotherapy units, use the following part numbers and descriptions.

	HENLEY No.	CGI No.
MODEL FLU110D & FLU110DE		
TOP SLEEVE	SLE0003	31775
TOP SLEEVE BEZEL	FRA0003	31456
ARM SLEEVE	SLE0002	31774
MODEL FLU115D & FLU115DE		
TOP SLEEVE	SLE0003	31775
TOP SLEEVE BEZEL	FRA0003	31456
INNER SLEEVE RETAINER	TRI0034	31883
SLEEVE RING	PPL0013	31679
OUTER SLEEVE RETAINER	TRI0035	31834
ARM SLEEVE	SLE0001	31773
CELLEX Dry Heat Medium 10 lbs (4.5 kg)		MED0001

MEDIUM

Use only Chattanooga's Cellex Dry Heat Medium in the Fluidotherapy unit. The Cellex medium is designed specifically for use in the Chattanooga Fluidotherapy units to ensure optimal and efficient operation of all the Fluidotherapy products.

WARRANTY

Fluidotherapy® - Dry Heat Therapy

DJO ("Company") warrants that the Fluidotherapy Dry Heat Therapy units ("Product") are free of defects in material and workmanship. This warranty shall remain in effect for one year (12 months) from the date of original consumer purchase. If this Product fails to function during the one year warranty period due to a defect in material or workmanship, at the Company's option, the Company or the selling dealer will repair or replace this Product without charge within a period of thirty (30) days from the date on which the Product is returned to the Company or the dealer.

All repairs to the Product must be performed by a service center certified by the Company. Any modifications or repairs performed by unauthorized centers or groups will void this warranty.

The warranty period for accessories is 90 days. Accessories consist of sleeves and use replaceable intake filter(s).

The warranty period for the motor(s) and distributor is 180 days.

To participate in warranty coverage, this Product's warranty registration card (included with Product) must be filled out and returned to the Company by the original owner within ten (10) business days of purchase.

This Warranty Does Not Cover:

Replacement parts or labor furnished by anyone other than the Company, the selling dealer, or a certified Company service technician.

Defects or damage caused by labor furnished by someone other than Company, the selling dealer, or a certified Company service technician.

Any malfunction or failure in the Product caused by product misuse, including, but not limited to, the failure to provide reasonable and preventive maintenance or any use that is inconsistent with the Product User's Manual.

COMPANY SHALL NOT BE LIABLE IN ANY EVENT FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES.

Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

To Obtain Service From Company or the selling dealer under this warranty:

1. A written claim must be made within the warranty period to the Company or the selling dealer. Written claims made to the Company should be sent to:

DJO, LLC
1430 Decision St.
Vista, CA 92081 USA
Phone: 1-800-592-7329 USA
Phone: 1-423-870-2281 or 1-317-406-2250
Fax: 1-317-406-2014

and

2. The Product must be returned to the Company or the selling dealer by the owner.

This warranty gives you specific legal rights and you may also have other rights which vary from state to state.

The Company does not authorize any person or representative to create for it any other obligation or liability in connection with the sale of the Product. Any representation or agreement not contained in the warranty shall be void and of no effect.

**THE FOREGOING WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED,
INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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