

Treadmill 4.0T

SPIRIT
powered by *dyaco*



The 4.0T treadmill is designed to provide optimal rehabilitation and fitness outcomes. Equipped with an extensive range of features, users will elevate their therapy sessions to new heights.

This treadmill offers unparalleled customization to cater to each client's specific needs:

- Full length handrails for added stability
- -5% - 15% decline/incline range
- Low starting speed for easy on/off access
- 20 innovative programs for a variety of patient's needs
- Current Leakage Protection with EMC Filter

The 4.0T is built to endure rigorous workouts and therapeutic sessions and is the ultimate solution for any commercial therapy environment.



Intuitive, easy-to-read display



Front and rear motors allow for true up/down walking



Direct speed and incline buttons for adjustments on the go

Specifications

| | |
|----------------------------|---|
| Incline | -5% to 15% |
| Speed | 0.3 to 12 mph |
| Motor | 4HP continuous duty |
| Belt and Deck | 22" x 60" self-lubricating deck |
| Step-Up Height | 7.75" |
| Power | AC 100-120V for North America AC 200-240V for EMEA |
| Max User Weight | 400 lb / 180 kg |
| Dimensions | 84.1" x 61.6" x 35.7" |
| Product Weight | 434.3 lb |
| Current Leakage Protection | Yes, with EMC Filter |
| Compliance | CE Class SA UL (Commercial grade) EN957 |
| Warranty | 5 Years parts, 1 Year labor |