

Upright Bike 7.0U

SPIRIT
powered by *dyaco*



Look no further for the best low impact lower body solution on the Rehabilitation market. Discover the revolutionary 7.0U upright bike, designed to redefine rehabilitation with unmatched comfort and precision.

Achieve your rehabilitation goals with exercises tailored to meet each client's specific needs:

- 50 levels of bidirectional resistance for a versatile workout experience
- A variety of preset programs, along with the option to create custom routines
- Fully adjustable high-density foam seat for maximum comfort
- Adjustable crank arms to accommodate any range of motion

The upright design promotes proper posture and alignment, enhancing the effectiveness of each training session.



Adjustable crank length to fit range of motion



Intuitive and easy-to-read display for a seamless workout experience



Numbered seat for easy adjustments

Specifications

Resistance System	Dual, electromagnetic resistance for constant and isokinetic modes
Drive Train	Heavy-duty, 8-groove poly V belt with spring-loaded idlers
Power	100-240 volts AC (standard power supply)
Workload Range	5 to 1700 Watts
Recommended User Height	4'10" - 6'6"
Max User Weight	440 lb / 200 kg
Dimensions	43.3" x 22.4" x 54.3"
Product Weight	117.5 lb
Compliance	CAN/CAS-C22.2 No. 60601-1:14 ANS/AAMI ES60601-1:2005+A2 (R2012) +A1 IEC 60601-1-2:2014 EN 60601-1-2:2015 IEC 60601-1:2015+A1:2012 EN 60601-1:2006+A1:2013+A12:2014 EN ISO 20957:2013 MDD 93/42/EEC Class Im, NB: 0123
Warranty	5 years parts, 1 year labor