

Recumbent Bike 7.0R

SPIRIT
powered by *dyaco*



Engineered to redefine rehabilitation with comfort and precision, the 7.0R recumbent bike sets the standard for advanced therapeutic solutions.

Experience unparalleled versatility with features that support clients at all stages:

- Swiveling seat for easy entry and exit
- 50 levels of resistance
- 11 preset programs for a variety of patient's needs
- Supportive high-density foam seat with fore/aft/recline angle adjustments

With these cutting-edge features, our 7.0R recumbent bike is your trusted partner in promoting rehabilitation, restoring mobility, and enhancing overall well-being.



Recumbent position reduces strain on joints and muscles



Contact heart rate monitoring capabilities



Forward and reverse pedaling motion

Specifications

Resistance System	Dual, electromagnetic resistance for constant and isokinetic modes
Drive Train	Heavy-duty, 8-groove poly V belt with spring-loaded idlers
Power	100-240 volts AC (standard power supply)
Workload Range	5 to 1700 Watts
Recommended User Height	4'10" - 6'6"
Max User Weight	440 lb / 200 kg
Dimensions	55.1" x 32.7" x 46.9"
Product Weight	180 lb
Compliance	CAN/CAS-C22.2 No. 60601-1:14, ANS/AAMI ES60601-1:2005+A2 (R2012) +A1 IEC 60601-1-2:2014, EN 60601-1-2:2015, IEC 60601-1:2015+A1:2012, EN 60601-1:2006+A1:2013+A12:2014, EN ISO 20957:2013, MDD 93/42/EEC Class Im, NB: 0123
Warranty	5 years parts, 1 year labor