

Recumbent Stepper 7.5S

SPIRIT
powered by *dyaco*



Make rehabilitation accessible for all with the 7.5S stepper, a semi-recumbent designed for wheelchair-bound patients. This stepper aligns with ADA compliance to facilitate accessibility and inclusivity in medical settings.

Patients with even the most limited abilities will benefit from features such as:

- Direct wheelchair access or 360-degree swiveling seat
- Bilateral movement and 20 levels of resistance
- 10-210 steps/min speed range for a variety of patient's needs
- 1"-12.5" adjustable step range

This innovative stepper offers low-impact, full-body exercise to aid in rehabilitation and improves cardiovascular health. Its adjustable settings ensure a safe and personalized experience for individuals of varying abilities.

**Optional MA902 Leg Stabilization Kit (Pair), Not Included*



Adjustable seat angle reduces back strain and ensures patient comfort



Full featured console with symmetry program, cooling fans and wireless heart rate



Direct wheelchair access

Specifications

| | |
|-------------------------|--|
| Braking Device | Eddy current brake with small flywheel for low starting inertia |
| Drive Train | Heavy-duty, 8-groove poly V belt with high-end linear bearings and guides |
| Power | 100-240 volts AC (standard power supply) |
| Workload Range | 5 to 1000 Watts |
| Recommended User Height | 4'10" - 6'6" |
| Max User Weight | 440 lb / 200 kg |
| Dimensions | 77.6" x 33.5" x 46.5" |
| Product Weight | 305.8 lb |
| Compliance | CAN/CAS-C22.2 No. 60601-1:14 ANS/AAMI ES60601-1:2005+A2 (R2012) +A1 IEC 60601-1-2:2014 EN 60601-1-2:2015 IEC 60601-1:2015+A1:2012 EN 60601-1:2006+A1:2013+A12:2014 EN ISO 20957:2013 MDD 93/42/EEC Class Im, NB: 0123 |
| Warranty | 5 years parts, 1 year labor |