


## DUAL COMFORT CORPAK™ WRIST WRAP HOT OR COLD THERAPY PACK

- Unique design comfortably contours to wrist and hand
- Strapping system provides compression and holds pack in place
- Soft and flexible, even when frozen
- Nontoxic, biodegradable fill maintains therapeutic temperature for 20 minutes
- Frost free exterior is gentle on the skin
- Plush surface slowly releases hot or cold therapy for subacute injuries
- Soft surface quickly releases hot or cold therapy for acute injuries



Available From:



### EASE PAIN AND DISCOMFORT

The Dual Comfort CorPak™ Wrist Wrap provides targeted hot or cold therapy as a versatile solution to ease pain and discomfort of the wrist and hand. Ideal for wrist and hand conditions such as carpal tunnel syndrome, arthritis, neuropathy, as well as acute or subacute injuries.

### COLD THERAPY

Use the Wrist Wrap cold to help reduce inflammation and swelling; providing soothing relief for muscle and joint pain. Conveniently store pack in the freezer for cold therapy the moment you need it! Help reduce inflammation caused by injuries such as strains, sprains, and other soft tissue injuries that may occur during activities. Also ideal for tendonitis, chronic pain, and bursitis.

### HOT THERAPY

Use the Wrist Wrap hot to accelerate recovery and healing by improving circulation and blood flow to the wrist and hand. Increasing the temperature of the afflicted area even slightly can soothe and relax muscles, promote healing of damaged tissue, and increase muscle flexibility.

### PROVIDES MILD COMPRESSION

The adjustable strap provides mild compression and secures the pack to the wrist and hand while keeping the fingers free. The pack's exterior is made of a proprietary blend of synthetic fibers and organic plant materials.

### SPECIFICATIONS

**Product Selection:** Dual Comfort CorPak™ Wrist Wrap Hot or Cold Compression Therapy Pack ACC-536

**Product Weight:** 2.19 lbs

**Dimensions:** Pack: 8.5" x 14" | Strap: 2" x 28.5"

**Materials:** Food grade propylene glycol USP, Bentonite clay, Water  
Cover: Wood pulp, Polyester, Nylon, Polyolefin alloy film  
Strap: Nylon hook, NEO-1/8 BKU/BKN

**Case Quantity:** 6

**Case Dimensions:** 12" x 12" x 6"

**Case Weight:** 15 lbs

**Latex:** This product is not made with natural rubber latex

**Country of Origin:** Made in USA

California Residents Only:

**⚠️WARNING:** Cancer and Reproductive Harm – www.P65Warnings.ca.gov