

hydraulic pinch gauges

 The Baseline® hydraulic pinch gauge uses a hydraulic system to assure measurement accuracy, repeatability and product reliability

- · Use for tip, key and palmar tests
- · Therapist can support the pinch gauge during testing, yielding a more accurate result for all pinch tests
- · Manufactured and serviced in the USA Lightweight units include protective carry case, instructions for use and norms
- Results consistent with published Baseline® and Jamar® studies
- · Maximum reading remains until reset
- Readings in lb and kg
- Lightweight unit includes protective carrying case
- CE certified

















| Model | LiTE® | Standard | HD® | Hi-Res [®] ER™ | Digital |
|------------------|------------------------|--------------------|------------------|---|---|
| Reading Capacity | 50 lb (23 kg) | 50 lb (23 kg) | 50 lb (23 kg) | 100 lb (45 kg) | 100 lb (45 kg) |
| Warranty | 1 year | 2 year | 5 year | 2 year | 2 year |
| Dial Size | 2.5" diameter | 2.5" diameter | 2.5" diameter | 3.5" diameter | 2.5" LCD screen |
| Color | Red | Blue | Black | Blue | Blue, black head |
| Notes | Inexpensive & reliable | Most popular model | Maximum warranty | Extra large dial, higher reading capacity | Toggle between max and current reading; Max clear button; 2 AAA batteries included |
| Sku | 12-0226 | 12-0235 | 12-0222 | 12-0228 | 12-0237 |

mechanical pinch gauges

- The Baseline® mechanical pinch gauge measures tip, key and palmar pinch strength in lb and kg
- Results are consistent with published Mathiowetz studies
- Indicator remains at the maximum reading until reset
- Includes protective case









| To the second |
|---------------|
| |

| Reading Capacity | 2 lb (0.9 kg) | 10 lb (4.5 kg) | 30 lb (13.6 kg) | 60 lb (27 kg) | |
|------------------|-------------------------------|----------------|--|---------------|--|
| Warranty | 2 year | 2 year | 2 year | 2 year | |
| Color | Gold | Silver | Blue | Red | |
| Notes | Use for weak or damaged hands | | Use for orthopedic and sports medicine | | |
| Sku | 12-0202 | 12-0203 | 12-0200 | 12-0201 | |



FABRICATION ENTERPRISES INC

PO Box 1500, White Plains, NY 10602 (USA) tel: 800-431-2830 / +1-914-345-9300 fax: 800-634-5370 / +1-914-345-9800 FabEnt.com info@fab-ent.com