




MMT manual muscle tester grip strength dynamometer

*Take manual muscle test readings
with your grip strength dynamometer*

Convert any
grip strength dynamometer

to a

MMT dynamometer



dynamometer, MMT
handle and 3 push pads

- Use the MMT handle to convert **ANY** hydraulic grip strength (hand) dynamometer to a MMT (Manual Muscle Test) dynamometer
- MMT handle comes with 3 push pads: straight, curved and small circular
- Use optional dual-grip handle to stabilize dynamometer during test procedure

MMT handle with 3 push pads

12-0580 fits all hydraulic dynamometers \$125.00

MMT / grip strength dynamometer w/ 3 pads

12-0590 LiTE 1 year warranty \$350.00

12-0591 standard 2 year warranty \$430.00

12-0592 HD 5 year warranty \$480.00

MMT / grip strength dynamometer w/ 3 pads and stabilizer handle


12-0593 LiTE 1 year warranty \$450.00

12-0594 standard 2 year warranty \$530.00

12-0595 HD 5 year warranty \$580.00

Dual-handled stabilizer

12-0583 fits all hydraulic dynamometers \$125.00



dynamometer, MMT
handle, stabilizer handle

