

# AUTO WAVE 695

## Radial Pressure Pulse Therapy



Mettler Electronics Corp.

# AUTO WAVE 695

## Radial Pressure Pulse Therapy



### Features

- Silent built-in compressor
- 7-color touch screen
- Ergonomic hand piece
- Treatment protocols with images
- Multiple language capability



### Indications

- Helps relieve minor muscle aches and pains
- Temporarily increases local blood circulation



### Applications for:

- Orthopedics
- Physical Therapy
- Podiatry
- Chiropractic







**Radial Pressure Pulse Therapy** is also known as Extracorporeal Pulse Activation Therapy. The Auto\*Wave 695 has a compressed air ballistic pulse energy generator that transfers pulse energy to the transmitter in the hand piece.

The motion and weight of the projectile, which is accelerated by compressed air, produces kinetic energy that is converted into pulse energy when the projectile strikes the surface of the transmitter.

The acoustic pulse energy is transformed into impact energy from the contact area into the therapeutic target tissue radially.

## ✿ Specifications

- **Pressure Range:** 1.0 - 4.0 Bar (0.1 Bar steps)
- **Frequency:** 1 - 17 Hz (1 Hz steps)
- **Pre-Set Protocols:** Yes
- **User-Defined Protocols:** 5,000
- **Dimensions:** 15.5" (L) x 15.3" (W) x 7.6" (H)
- **Weight:** 18.7 lbs (Headpiece excluded)



Ballistic Mechanism of Radial Pressure Pulse Therapy

## ✿ Transmitter

Transmitters						
Type	6 A (D = 6mm Acupuncture)*	15 R (D = 15mm Radial)	15 F (D = 15mm Focus)*	15 D (D = 15mm Deep)*	20 R (D = 20mm Radial)	36 R (D = 36mm Radial)
Penetration Depth	0 - 40 mm	0 - 40 mm	0 - 20 mm	0 - 60 mm	0 - 50 mm	0 - 50 mm

\* Optional Accessories

## Certification

CE 0537

CE (MDD)



Health Canada



### New Custom Cart

- Four locking wheels
- Large handle to easily move the cart
- Three internal shelves to store accessories and supplies
- Easy opening magnetically secured door



**Mettler Electronics Corp.**

**Mettler Electronics Corp.**

**1333 South Claudina Street**

**Anaheim, CA 92805 U.S.A.**

**Toll Free: 1.800.854.9305**

**Tel: 1.714.533.2221**

**Fax: 1.714.635.7539**

**Email: [mail@mettlerelectronics.com](mailto:mail@mettlerelectronics.com)**

**[www.mettlerelectronics.com](http://www.mettlerelectronics.com)**