



The RecoveryAir JetBoots Solution

RecoveryAir JetBoots

Integrated Pneumatic Compression System

Fully wireless compression boots.

Portable. Convenient. Fast.

Pneumatic compression is more convenient than ever while still allowing you to experience the full benefits of RecoveryAir’s advanced technologies. It eliminates the external pump and tubes, replacing them with integrated motors, making the JetBoots ultra-portable. The RecoveryAir JetBoots’ components are all packed into a sleek, compact design for easy storage and use anywhere.

Recovery Jetboot Difference

TruGrade Technology: RecoveryAir’s TruGrade Technology™ creates precise pressure, starting from your feet moving toward your heart, to temporarily increase circulation to the treated area.

FastFlush Technology: Aids in recovery by completing one full cycle in just 60 seconds. 2-3x faster than competitors

What Is Pneumatic Compression?

RecoveryAir’s pneumatic compression is “mechanical massage”, built to boost circulation, decrease swelling, and reduce muscle soreness and fatigue.

- Benefits:**
- Boosts circulation
 - Reduces muscle soreness
 - Decreases swelling and stiffness
 - Relieves muscle fatigue

Features:

60 second full flush
FastFlush Technology



4 internal overlapping chambers
Superior, hygienic design



Precise pressure control
Safe & effective



240 minutes of battery life
Convenient & portable



2 customizable programs
Self-guided treatments



Additional Features:

- 4 time intervals (20, 40, 60, and continuous)
- Basic Program Sharing
- My Recovery Stats Tracking
- Bluetooth connectivity

Specs:

Dimensions (per single boot):

- Small - 36.8 in (L) x 5.2 in (W) x 13.7 in (H)
- Medium - 43.3 in (L) x 5.2 in (W) x 15.1 in (H)
- Large - 46.5 in (L) x 5.2 in (W) x 15.1 in (H)

Weight (per single boot):

- Small - 4.12 lbs
- Medium - 4.40 lbs
- Large - 4.48 lbs

Precision:

Adjustable Pressure Range (25-100 mmHg) in increments of 25 mmHg*

Gradient Type:

Full Negative

Decompression:

Full and Rapid Release

Cycle Time:

60 second full flush (2-3x the speed of competitors)

Material:

Medical-grade coated fabric

*Adjustable pressure changes via the Therabody app in increments of 5 mmHg

What’s Included:

- 2 compression boots with integrated pumps
- 36W Charger - Jack Port
- Splitter Charging Cord
- Boots Pouch

Compression Boots Sizing Guide

Boot Size	Height	Length Range	Max Width
Small	5’-5’5” 152-169cm	27”-33.5” 70-85cm	27.6” 70cm
Medium	5’6”-6’ 170-183cm	31.5”-37.5” 80-95cm	31.5” 80cm
Large	6’-6’6” 183-198cm	35.5”-41” 90-105cm	31.9” 81cm

If you are in-between sizes, we recommend sizing down to ensure the best fit