

WORLD'S LIGHTEST HEAVY DUTY WHEELCHAIR!

# Feather CHAIR HD

Where heavy duty meets light weight.



The Feather Chair HD boasts a 22" wide seat, 350 lb. capacity and weighs only 22 lbs., **Pop off the wheels**, and it weighs only 15 lbs.! Swing-away footrests and anti-tippers offer comfort and safety. The **back folds down** for extra added convenience during storage or travel.

## FEATURES:

- only 15 lbs without wheels
- 22" wide seat
- folding frame
- anti-tippers
- optional elevating leg rests

## SPECIFICATIONS

Weight Capacity	350 lbs.		
Product Weight	22 lbs.	Front Wheel Size	8"
Weight Without Wheels	15 lbs.	Rear Wheel Size	24"
Seat Width	22"	Axle Type	Fixed, Quick Release
Seat Depth	16"	Armrest Options	Full Length
Seat-to-Floor Height	19"	Armrest Height	9"
Overall Width	32"	Back Height	16"
Overall Length	34"	Upholstery Type	Nylon
Overall Height	28"	Upholstery Color	Black w/ Green Trim
Folded Depth	13"	Footrests	Swingaway
Folded Height	28"	Warranty on Frame	Lifetime Limited
Folded Width	15"	Warranty on Wearables	Six Months Limited

Lightweight



Foldable



Locking Brakes



Arrives Assembled

