



Bariatric Drop-Arm Commode

The **ProBasics Bariatric Drop-Arm Commode** features a durable steel frame with an easy-release push button arm mechanism allowing for safe lateral patient transfers to and from the commode. Equipped with wider seat and reinforced steel frame to support patient weights up to 650 pounds, it comes complete with commode bucket and lid.

FEATURES

- Features a reinforced steel frame and an extra wide platform seat designed for comfort and stability
- Equipped with an easy-to-release arm mechanism allowing for safe lateral patient transfers to and from the commode
- Bucket easily slides into place
- Contoured armrests provide additional patient comfort
- Comes complete with a commode bucket and cover
- Color: gray frame, gray accessories
- HCPCS Code: E0168*

SPECIFICATIONS

Width between arms:	24.5"	
Outside width:	28.5"	
Seat dimensions:	23"W x 18"D	
Seat opening:	8.25"W x 10"D	
Adjustable seat height:	19.5" – 23.5"	
Weight Capacity:	650 lb	
Shipping Dimensions:	29.5" x 21" x 20"	
Shipping Weight:	26 lb	
Warranty:	Limited lifetime on frame One year on all parts	

Item #	Description	UPC	UOM
BSBDAC	ProBasics Bariatric Drop-Arm Commode with 24" Platform Seat	815067071869	1/cs

^{*} Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.