

CANDO® Web

Finger Flexion, Extension & Wrist Exerciser

Tear-resistant formula

- Perform flexion, extension, opposition and supination exercises
- Modify exercises by adjusting the hand position or depth of finger insertion
- Modify resistance by switching to a different color-coded resistance web
- Available in both 14" and 7" diameter size; latex and latex-free
- Resistance levels are color-coded
- Multi-resistance latex webs have two resistance levels in one web



Perform flexion, extension and wrist exercises




Large 14" diameter

Small 7" diameter

		Large 14" diameter		Small 7" diameter	
		latex	latex-free	latex	latex-free
	xx-light	10-0850	10-0870	10-0860	10-0880
	x-light	10-0851	10-0871	10-0861	10-0881
	light	10-0852	10-0872	10-0862	10-0882
	medium	10-0853	10-0873	10-0863	10-0883
	heavy	10-0854	10-0874	10-0864	10-0884
	x-heavy	10-0855	10-0875	10-0865	10-0885
	set of 6	10-0856	10-0876	10-0866	10-0886

Multi-resistance 14" diameter, latex

		xx-light / medium	10-0857
		light / heavy	10-0859



Multi-resistance 14" large diameter, latex

Compare then decide

CanDo® Web

Power Web®

Resistance levels available 14"

6 standard, 2 multi level, 6 latex-free

6 standard, 4 multi level, 1 latex-free

Resistance levels available 7"

6 standard, 6 latex-free

3 standard, 1 latex-free

Material

Newly formulated to be more tear resistant

Design

Round finger holes

Square finger holes



FABRICATION ENTERPRISES INC
 PO Box 1500, White Plains, NY 10602 (USA)
 tel: 800-431-2830 / +1-914-345-9300
 fax: 800-634-5370 / +1-914-345-9800
 FabEnt.com info@fab-ent.com