

# Clinical Monofilaments

- Same trusted product in new color-coded system
- Measure and detect skin sensitivity, nerve damage or impairment, and changes in neurological status
- Color code provides an easy reference for recording results: green, normal; blue, diminished light touch; purple, diminished protective sensation; red, loss of protective sensation
- **Fold-Up™ Monofilaments** have a top that folds to protect monofilament without an additional case
- **Tactile™ Monofilaments** come with protective case

### Size (grams of force)

- 1.65 (0.008 gm), 1 each
- 2.36 (0.02 gm), 1 each
- 2.44 (0.04 gm), 1 each
- 2.83 (0.07 gm), 1 each
- 3.22 (0.16 gm), 1 each
- 3.61 (0.4 gm), 1 each
- 3.84 (0.6 gm), 1 each
- 4.08 (1.0 gm), 1 each
- 4.17 (1.4 gm), 1 each
- 4.31 (2.0 gm), 1 each
- 4.56 (4.0 gm), 1 each
- 4.74 (6.0 gm), 1 each
- 4.93 (8.0 gm), 1 each
- 5.07 (10.0 gm), 1 each
- 5.07 (10.0 gm), 100 each
- 5.18 (15.0 gm), 1 each
- 5.46 (26.0 gm), 1 each
- 5.88 (60.0 gm), 1 each
- 6.10 (100.0 gm), 1 each
- 6.45 (180.0 gm), 1 each
- 6.65 (300.0 gm), 1 each

color
normal
diminished light touch
diminished protective sensation
loss of protective sensation



### Fold-Up™

- 12-1720
- 12-1721
- 12-1722
- 12-1723
- 12-1724
- 12-1725
- 12-1726
- 12-1727
- 12-1728
- 12-1729
- 12-1730
- 12-1731
- 12-1732
- 12-1733
- 12-1733-100
- 12-1734
- 12-1735
- 12-1736
- 12-1737
- 12-1738
- 12-1739

### Tactile™

- 12-1640
- 12-1641
- 12-1642
- 12-1643
- 12-1644
- 12-1645
- 12-1646
- 12-1647
- 12-1648
- 12-1649
- 12-1650
- 12-1651
- 12-1652
- 12-1653
- 12-1653-100
- 12-1654
- 12-1655
- 12-1656
- 12-1657
- 12-1658
- 12-1659

### Clinical Monofilament Sets

- 2-piece, normal, 2.83, 3.61
- 2-piece, protective sensation 4.56, 5.07
- 5-piece, hand, 2.83, 3.61, 4.31, 4.56, 6.65
- 6-piece, foot, 2.83, 3.61, 4.31, 4.56, 5.07, 6.65
- 20-piece, complete

- 12-1745
- 12-1746
- 12-1740
- 12-1741
- 12-1743

- 12-1660
- 12-1668
- 12-1662
- 12-1664
- 12-1666

### Monofilament Accessories

- Use to keep track of evaluation results; 1 pad with 100 sheets
- 12-1360 hand sheet, 1 pad
- 12-1361 foot sheet, 1 pad
- 12-1362 5 pencils (R, G, B, P, BLK)

### Baseline® Fold-Up™ Monofilaments



20-piece Fold-Up™ set complete



### Baseline® Tactile™ Monofilaments



6-piece Tactile™ set, foot

