

ComforTrac Cervical Traction



The ComforTrac® cervical home traction unit has been carefully engineered to provide safe and effective traction to the cervical region.

The design ensures maximum comfort, efficacy, and ease of use. The ComforTrac is sturdy enough to be used in a clinical setting or the home. Therapists can either sell or rent units.

Key Features:

- Up to 23 kg (50 lbs) of cervical traction
- 3 Angles of incline: 10, 15 and 20 degrees (with no additional parts or attachments)
- Memory pillow conforms comfortably to the head
- Independently adjustable neck wedges provide a custom fit
- Hand pump makes it easy to increase or decrease traction force; Displays force in both kilograms and pounds
- Quick release button releases traction force
- Lightweight durable carry-case
- 3-year warranty

50-1070 ComforTrac home cervical traction

SUPPORTING RESEARCH AND REFERENCES:

Spine. 1996 Aug 15;21(16):1877-83. "Nonoperative management of herniated cervical intervertebral disc with radiculopathy." PMID: 8875719

Spine. 1985 Dec;10(10):867-71. "Cervical spine disorders. A comparison of three types of traction." PMID: 3914085

J Manipulative Physiol Ther. 2002 Mar-Apr;25(3):188-92. "Intermittent cervical traction for cervical radiculopathy caused by large-volume herniated disks." PMID: 11986581

Bangladesh Med Res Counc Bull. 2002 Aug;28(2):61-9. "Effects of cervical traction and exercise therapy in cervical spondylosis." PMID: 12825763



ComforTrac Lumbar Traction



The ComforTrac® lumbar home traction unit is carefully engineered to provide maximum comfort, efficacy, and ease of use.

Key Features:

- Replicates clinical traction
- Delivers up to 90 kg (200 lbs) of lumbar traction
- Friction-free traction provides smooth, consistent and even force
- Hand pump makes it easy to increase or decrease traction force; Displays force in both kilograms and pounds
- Quick release hand pump button releases traction force
- Automatic carriage return
- Padded surface provides superior comfort
- Non-slip harness belt and durable belt clips secure patient during traction
- Transporting the unit is effortless; weighs 6.5 kg (14 lbs)
- Carry case with handle is included
- No assembly required

Specifications:

- Product weight: 6.5 kg (14 lbs)
- Shipping weight: 8.5 kg (19 lbs)
- Dimensions: 63.5 cm x 89 cm x 7.6 cm (25" x 35" x 3")

50-1075 ComforTrac home lumbar traction

SUPPORTING RESEARCH AND REFERENCES:

Komori H, Shinomiya K, Nakai O, et. al., The Natural History of Herniated Nucleus Pulposus with Radiculopathy. Spine 21:225-229, 1996.

Onel D, Tuzlaci M, Sari H, et. al. Computed Tomographic Investigation of the Effect of Traction on Lumbar Disc Herniations. Spine 14:82-90, 1989.

Matthews JA, Hickling J, Lumbar Traction: A Double-Blind Controlled Study for Sciatica. Rheumatology and Rehabilitation 14:222-225, 1975.



pump and belt clips are patented



carry case