

Digi-Flex® Hand Exercisers

Select the hand exerciser that's best for you!

- Use to rehabilitate injuries or conditions affecting the hand and upper extremities
- Develops isolated finger strength, flexibility, and coordination
- Available in eight (8) strength levels
- Comes with exercise pamphlet

Digi-Flex®

- Digi-Flex® model has upper and lower spring to allow for more hand motion flexibility. Use for individual finger or hand



Digi-Flex LiTE®

- Digi-Flex LiTE® is An inexpensive option for patient giveaway or sale. It has individual finger locks and a single spring design



patented

Digi-Flex Thumb®

- Digi-Flex Thumb® allows for thumb involvement. The thumb piece can attach at three different positions to accommodate all thumb sizes and swivels for all thumb positions



patented

force to compress
finger/hand (lb)

Color	Strength Level	Force to Compress (lb)	Weight (lb)	Weight (kg)	Product Code
tan	xx-light	0.75	2.5		10-0751
yellow	x-light	1.5	5.0		10-0740
red	light	3.0	10.0		10-0741
green	medium	5.0	16.0		10-0742
blue	heavy	7.0	23.0		10-0743
black	x-heavy	9.0	30.0		10-0744
silver	xx-heavy	11.0	36.0		10-0752
gold	xxx-heavy	13.0	43.0		10-0753

yel	red	grn	blu	blk	set of 5 (1 each)	10-0745			
yel	red	grn	blu	blk	set with plastic stand	10-0746			
yel	red	grn	blu	blk	set with metal stand	10-0749			
tan	yel	red	grn	blu	blk	sil	gld	set of 8 (1 each)	10-0758
tan	yel	red	grn	blu	blk	sil	gld	set of 8 with 2 metal stands	10-0759

Digi-Flex®


Digi-Flex LiTE®

Digi-Flex Thumb®

10-3770
10-3771
10-3772
10-3773
10-3774
10-3775
10-3776
10-3777

10-3760
10-3761
10-3762
10-3763
10-3764
10-3765
10-3766
10-3767

10-3769
10-3785
10-3786
10-3768
10-3787



metal stand

Digi-Flex® Stands

- Holds up to 5 Digi-Flex® units
- | | |
|---------------|---------|
| metal stand | 10-0750 |
| plastic stand | 10-0747 |



FABRICATION ENTERPRISES INC
 PO Box 1500, White Plains, NY 10602 (USA)
 tel: 800-431-2830 / +1-914-345-9300
 fax: 800-634-5370 / +1-914-345-9800
 FabEnt.com info@fab-ent.com

Digi-Flex Multi® Hand Exerciser

Create a customized hand exerciser in 3 easy steps!



pinch release tabs and pull up to remove button



insert button until a "snap" is heard



use your customized hand exerciser

Clinic Sets

Customize hand therapy like never before. Adjust each finger's resistance to fit the patient's rehab program.

Clinic sets are available with either 5 bases or 5 pre-built exercisers, 20 or 32 piece resistive button sets with storage case, and rack.

More bases, exercisers, buttons, and button sets are available for purchase as needed. See the following link for more information: Fab-Ent.com/DigiFlexMulti

small clinic set with 5 bases


10-3840 deluxe, 32 buttons 4 ea tan yel red grn blu blk sil gld

10-3841 standard, 20 buttons 4 ea yel red grn blu blk

large clinic set with 5 pre-built exercisers

10-3845 deluxe, 32 buttons 4 ea tan yel red grn blu blk sil gld

10-3846 standard, 20 buttons 4 ea yel red grn blu blk



small clinic set with 5 bases, deluxe (10-3840)

patented



large clinic set with 5 pre-built exercisers, standard (10-3846)

Individual Finger Locks

Insert the white finger lock button when finger movement or exertion is counter-indicated.



10-3738-4 4 each with case



FABRICATION ENTERPRISES INC
 PO Box 1500, White Plains, NY 10602 (USA)
 tel: 800-431-2830 / +1-914-345-9300
 fax: 800-634-5370 / +1-914-345-9800
 FabEnt.com info@fab-ent.com