


BASELINE[®] EVALUATION INSTRUMENTS

REF 12-0297

PNEUMATIC BULB HAND DYNAMOMETER / PINCH GAUGE



Adjustable dynamometer/pinch gauge measures both grip and pinch strength. Comes with three different sized interchangeable squeeze bulbs that can be easily fitted to the gauge. Lightweight and portable. Maximum reading remains until reset by turning center knob. Unit is calibrated in PSI (0-30) and kPa (0-200).

⚠ This Product Contains Dry Natural Rubber.



Fabrication Enterprises Inc.
250 Clearbrook Road, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
FabEnt.com



Authorized CE representative:
AJW Technology Consulting GmbH
Königsallee 106
40215 Düsseldorf (Germany)