

Wall Growth Chart

User Manual

Intended Purpose

The Baseline[®] Wall Growth Chart is a stick-on measuring tool designed for parents and pediatricians to accurately track height measurements up to 7 feet for people of all ages.

Introduction

The Baseline[®] Wall Growth Chart is your go-to tool for tracking height quickly and accurately. Whether you're a parent documenting your child's growth or a pediatrician monitoring patient progress, this chart is designed to meet your needs effortlessly. It precisely measures height up to 7 feet (0 to 84 inches). It offers a clean and simple design that fits seamlessly into any space.

Extremely portable and easy to install, simply place it along any wall starting at floor level for instant use. Add it to a wall at the doctor's office or at home to keep a record of your child's growth milestones. The blue text ensures easy readability, even in low-light settings. The Baseline[®] Wall Growth Chart is a reliable tool for patient height tracking and is suitable for all ages.

How to Install

The Baseline[®] Wall Growth Chart has an adhesive back. Peel and stick the wall growth chart to your desired wall. Start from the floor, or cut to your desired starting point.

Measurement Tolerance

This device measures inches and millimeters with a tolerance of ±0.04 inches (±1 millimeter)

Disposal Method

Dispose of item in accordance with the local/regional/ national/international regulations.

⚠ WARNING:

- Verify installation height and alignment
- Maintain proper lighting for accurate readings
- Replace if measurement markings become unclear
- Keep away from direct sunlight to prevent fading
- Ensure stable wall mounting
- Verify accuracy of reference points

