

**BASELINE[®] EVALUATION
INSTRUMENTS**
**plastic 6-inch 180° pocket
goniometer**

REF
12-1005



Intended Purpose: permits observation of joint's axis of motion, and range of motion

- 180° head has three scales for ease of use
- primary scale reads in 1° increments
- goniometer arm has a linear scale in inches and centimeters
- goniometer arm remains at measurement until reset

other goniometer styles and sizes available:

12" 360°		12-1000*
8" 360°		12-1001*
6" 360°		12-1002*
6" 180° pocket		12-1005*
6" 360° rulometer	...	12-1006*
3" -30° to 120° finger	...	12-1012*

Measurement Tolerance
This device measures inches and centimeters with a tolerance of ±0.04 inches (±1 millimeter) and degrees with a tolerance of ±1 degree.

Disposal Method
Dispose of item in accordance with the local/ regional/national/international regulations.

WARNING:

- Clean and disinfect device before and after each patient use
- Verify measurement markings are clearly visible
- Ensure smooth operation of all moving parts
- Align device properly with anatomical landmarks
- Apply gentle pressure during measurements
- Stop if patient reports pain or discomfort
- Store flat to prevent warping
- Replace if accuracy becomes compromised
- Do not attempt to bend or modify the device


FEI
Fabrication Enterprises Inc.
250 Clearbrook Road, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
FabiEnt.com


CE
2862


EC REP
AJW Technology Consulting Center
Breits Strasse 3
40213 Düsseldorf Germany


MD


Mediway
HL-18-000000248
Prinses Margrietplantsoen 33, Suite 123
2055 AM, The Hague (NL)

©2025 FEI, all rights reserved.
Baseline and Hi-Res are
trademarks of Goldberg. rev0325