



plastic digit goniometer



LOT

12-1014



Readings can be taken by the therapist using only one hand; this frees up the therapist's other hand to stabilize the patient during measurement or to record the measurement.

Measures 110° flexion through 40° hyper-extension in 2° increments.

Arm has a 10cm/4" scale.

Made of sturdy ABS plastic with white markings on black background.



Fabrication Enterprises Inc.
250 Clearbrook Rd, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
FabEnt.com

EC | **REP**



AJW Technology Consulting GmbH
Königsallee 106
40215 Düsseldorf (Germany)