


# Feeder Seat Measurements



Size	REF	Inside seat width	Depth	Height from bottom of the seat to harness hole*
Small	30-1070	8.0" (20.32 cm)	6" (15.24 cm)	11" - 13" (27.94 - 33.02 cm)
Medium	30-1071	10.0" (25.4 cm)	8" (20.32 cm)	14" - 17" (35.56 - 43.18 cm)
Large	30-1072	12.5" (31.75 cm)	12" (30.48 cm)	19" - 23" (48.26 - 58.42 cm)
X-Large	30-1074	15.5" (39.37 cm)	14" (35.56 cm)	20" - 24" (50.8 - 60.96 cm)

\*Range accounts for the harness hole length. Measurements are taken from the bottom of the seat to the bottom of the harness hole and the top of the harness hole.