

# Retention Straps for Swing Seats


## User Guide

### Introduction

The Retention Straps for Swing Seats are designed to make entry and exit for swing seats easier for both the user and their caregiver. The straps anchor the seat to the frame, holding the seat steady to prevent swaying during entry and exit.

### Setup

- Attach retention strap end without clip to frame.
- Secure the other end of the retention strap to swing seat using clip.
- Adjust the length of the retention strap using the buckles.



Fabrication Enterprises Inc.  
250 Clearbrook Road, Suite 240  
Elmsford, NY 10523 (USA)  
tel: +1-914-345-9300 • 800-431-2830  
fax: +1-914-345-9800 • 800-634-5370  
FabEnt.com



Authorized CE representative:  
AJW Technology Consulting GmbH  
Königsallee 106  
40215 Düsseldorf (Germany)