

#### **Lo-Boy Whirlpools**



### Division of Acorn Engineering Company®



## L Series 90 Gallon Mobile Whirlpool

L-90-M



L-90-M

Whitehall L-90-M Mobile Whirlpools are recommended for treatment of extremities, hip and back area and allows space for exercising.

Tank Capacity is 90 gallons.

Whirlpool Tank is 60"L x 24"W x 18"D and is seam-less welded construction, fabricated from heavy gauge, type 304 stainless steel, polished to a satin finish. Tank rim is reinforced with a concealed stainless steel rod. Tank bottom is a seamless, coved (rounded) design for effective cleaning to minimize bacteria buildup.

**Heavy Duty Swivel Casters** support all mobile models. Rear pair of casters are lockable for stabilizing whirlpool during treatment.

**Turbine Raising and Lowering Device** permits adjustment of desired height and direction of water agitation towards the area of the body to be treated. It is spring balanced and provided with a locking device.

**Turbine Assembly** functions both as agitator and emptying device. It is equipped with Underwriters Laboratory listed 1/2 HP jet pump motor protected with an ABS plastic shell, automatic thermal overload protector, and lifetime-sealed bearings. All parts contacting the water are stainless steel, chrome plated brass or bronze.

**Drain option -SDP** has a bottom drain, leading to a separate, concealed water evacuation pump. A hose is attached to the side of the unit to facilitate the evacuation.

**Thermometer** is provided with a dual scale, adjustable calibration, and a stem inserts into a retainer to prevent rattle during whirlpool agitation.

#### **GUIDE SPECIFICATIONS**

Provide Whitehall Mobile Whirlpool (specify model number and options) with a 90-gallon capacity, 60"L x 24"W x 88"D. Whirlpool shall be fabricated from heavy gauge, type 304 stainless steel. Construction shall be seamless welded and exposed surfaces shall be polished to satin finish. Bottom of tank shall be a seamless, rounded coved design to minimize bacteria build up. Whirlpool shall have heavy duty swivel casters with locking device on rear pair of casters. Provide a turbine assembly with raising and lowering device which functions both as agitator and emptying device, permits adjustment of desired height, direction of water agitation, and can be locked into place. Turbine assembly shall be UL listed 1/2 HP jet pump motor with automatic thermal overload protector and lifetime-sealed bearings. Provide a thermometer with a stem retainer to prevent rattling when whirlpool is in operation.

Please visit www.whitehallmfg.com for most current specifications.





Whitehall Mfg. • P.O. Box 3527 • City of Industry, CA 91744 • (800) 782-7706 • (626) 968-6681 • Fax (626) 855-4862

Updated: 10/03/18



# **Lo-Boy Whirlpools**

MODEL NUMBER AND OPTION SELECTIONS	Suffix Options (See Whirlpool Options Section for Details)
Series ☑ L Lo-Boy Series	<ul><li>□ -ADH Auxiliary Drain Hose Assembly</li><li>□ -ALC Adjustable Low Chair</li><li>□ -ALS Adjustable Low Stool</li></ul>
Capacity ☑ -90 Gallons	☐ -ARS Arm Rest Support ☐ -AS3 Adjustable Suspension Seat
Type ☑ -M Mobile Whirlpool	☐ -EG Enviro Glaze Color Powder Coating - Specify Color:
Turbine Electrical Requirements ☐ 115VAC / 50Hz ☐ 115VAC / 60HZ ☐ 230VAC / 50Hz ☐ 230VAC / 60HZ	(Special Finishes Consult Factory)  ☐ -GFI Ground Fault Interrupter ☐ -HR Head Rest, Adjustable ☐ -PCS Pump Strainer Cover ☐ -PHS2 Patient Handling System (Chair Type)
Country of Use	-RHG Rubber Hand Grips
For Mixing Valve: Refer to supplemental MXT15 sheet for details.	☐ -RS Reclining Seat ☐ -SDP Separate Drain Pump ☐ -TTS2 Tank Top Seat
Amps: 6 9/6 2 Hertz: 50/60 Amps: 3 5/3 1 Hertz: 50/60 Medical	d per US & Canadian Standards; UL & CSA C22.2 No.  Certified per European Standards; CE CSA C22.2 No. 60601-1  Cortified per European Standards; CE CSA C22.2 No. 60601-1
	60" (1525) 5 (610) 6 4
MOBILE L-90-M S NOTE: 1. 1/2 hp Jet Pump Motor 2. Variable Pressure Control 3. Thermometer	SERIES WHIRLPOOL  4. Stainless Steel Tank 5. Auxiliary Drain Valve 6. Turbine Raising & Lowering Device
	WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov
All dimensions are nominal and subject to	& APPROVAL FOR MANUFACTURING
Whitehall assumes no responsibility fo use of void or superseded data © Whitehall Model Number & Options	_Quantity
International.	
Please visit <b>www.whitehallmfg.com</b> for	

Whitehall Mfg. • P.O. Box 3527 • City of Industry, CA 91744 • (800) 782-7706 • (626) 968-6681 • Fax (626) 855-4862

Approval for Manufacturing/Signature\_