



PROTECT YOURSELF FROM EVERYDAY GERMS AND BACTERIA

For more details visit:
www.theshield.pro



INSTALLATION INSTRUCTIONS

Step 1. Place TheShield.Pro over top of wheel and gently press down to fit. The logo should face out with the V shape toward the front of wheel. (Diagram #1)

Step 2. Attach front strap to the front of wheel chair frame using Velcro.

(Diagram #2)

Step 3. Thread back straps through arm rest, then attach Velcro behind chair.

(Diagram #3)



To propel yourself forward gently press down and push. TheShield.Pro should glide easily along with wheel chair propulsion. If TheShield.Pro grips too tightly, remove, review directions and reattach.

Manufactured by The Shield, LLC | Hatley, Wisconsin 54440

