

# cold therapy products

## English-style ice cap



- bags are made of waterproof, rubberized fabric with plastic screw caps

| diameter | each         | bulk (12 each)    |
|----------|--------------|-------------------|
| 6"       | 11-1060 5.00 | 11-1060-12 55.00  |
| 9"       | 11-1061 6.00 | 11-1061-12 67.50  |
| 11"      | 11-1062 9.00 | 11-1062-12 102.50 |

## Relief Pak® insulated ice bags



- refillable and disposable
- available with hook/loop band or string ties
- easy-to-fill neck opening
- small bag: 5" x 13"
- large bag: 7" x 13"

|       | hook/loop band | each | bulk (10 each)   |
|-------|----------------|------|------------------|
| small | 11-1241        | 4.00 | 11-1241-10 37.50 |
| large | 11-1240        | 4.00 | 11-1240-10 37.50 |
|       | tie strings    |      |                  |
| small | 11-1243        | 4.00 | 11-1243-10 37.50 |
| large | 11-1242        | 4.00 | 11-1242-10 37.50 |

## instant cold packs



- squeeze the bag and in just 3 seconds the pack becomes 34°F
- cold is achieved by safe chemical reaction
- disposable, single use

|         |         |       |       |
|---------|---------|-------|-------|
| 11-1020 | 6" x 9" | 12 ea | 15.00 |
| 11-1021 | 4" x 6" | 12 ea | 12.50 |

## AirCast® CryoCuff® cold compression

- combines compression and cold: compression minimizes swelling while cold minimizes pain
- cuffs made to fit specific body parts
- cooler holds enough ice and water for 8 hours of cryotherapy
- cuffs are made of nylon and vinyl
- cuffs should be used with gravity feed cooler or cooler with continuous AutoChill® electric pump (sold separately)



ankle



knee

### CryoCuff® (cuff and cooler)

| cuff and cooler |        |              | cuff only |        |
|-----------------|--------|--------------|-----------|--------|
| 11-1550         | 147.50 | ankle        | 11-1568   | 82.50  |
| 11-1552         | 147.50 | medium foot  | 11-1570   | 82.50  |
| 11-1556         | 147.50 | medium knee  | 11-1573   | 82.50  |
| 11-1557         | 147.50 | large knee   | 11-1574   | 82.50  |
| 11-1560         | 165.00 | shoulder     | 11-1577   | 100.00 |
| 11-1561         | 165.00 | XL shoulder  | 11-1578   | 100.00 |
| 11-1563         | 147.50 | XL thigh     | 11-1582   | 82.50  |
| 11-1564         | 147.50 | calf         | 11-1583   | 82.50  |
| 11-1565         | 147.50 | back/hip/rib | 11-1584   | 82.50  |
| 11-1566         | 147.50 | elbow        | 11-1585   | 82.50  |
| 11-1562         | 147.50 | thigh        | 11-1581   | 82.50  |
| 11-1567         | 147.50 | hand/wrist   | 11-1586   | 82.50  |

### accessories

|         |                                   |        |
|---------|-----------------------------------|--------|
| 11-1588 | cooler with tubing system         | 65.00  |
| 11-1589 | AutoChill® pump                   | 207.50 |
| 11-1549 | AutoChill® pump and cooler system | 205.00 |
| 11-1548 | IC cooler pump and cooler system  | 155.00 |



shoulder



thigh



back/hip/rib



wrist



cooler



AutoChill® pump