

Relief Pak® Cold n' Hot® Sensaflex® compress



- Cold n' Hot® Sensaflex® compress has a comfortable soft fabric that conforms to the curves of your body, is gentle on the skin and washable
- chill in freezer or heat in microwave
- when chilled, soothes discomfort from injury, sprains and strains, muscle soreness, bruising, swelling and tension headaches
- when heated, soothes discomfort from stiff neck, backache, arthritis, sprains and bruises
- non-toxic, non-crystallizing gel allows pack to contour to body
- standard hook and loop adheres to pack to hold pack in place



			each		bulk (case qty)	
F	circular	10" diameter	11-1297	12.50	11-1297-12 (12 ea)	142.50
G	small	3" x 5"	11-1291	7.50	11-1291-12 (12 ea)	82.50
H	medium	7" x 12"	11-1292	10.00	11-1292-12 (12 ea)	112.50
I	large	10" x 13"	11-1293	12.50	11-1293-6 (6 ea)	70.00
J	neck	7" x 22"	11-1294	12.50	11-1294-12 (12 ea)	142.50
K	tri-sectional	8" x 16"	11-1298	12.50	11-1298-6 (6 ea)	70.00



F



G



H



I



J




K



Relief Pak® Cold n' Hot® elastomer pack

- designed for moist heat and soothing cold applications
- if punctured, pack won't leak and will remain usable
- tough, flexible gel, specially selected fabric and velcro straps provide an exceptionally durable product
- pack maintains good flexibility at -20°F
- stretch cover allows conformity, heat transfer and comfort
- heat pack in microwave... NOT in water

pad		each		bulk (case of 12)	
L	12½" x 12½"	11-1631	37.50	11-1631-12	427.50
wraps					
M	5" x 9½"	11-1633	20.00	11-1633-12	227.50
N	7¼" x 17¼"	11-1634	32.50	11-1634-12	370.00
O	10" x 24"	11-1635	57.50	11-1635-12	655.00
neck/shoulder					
P	11½" x 20"	11-1636	42.50	11-1636-12	482.50
head					
Q	3" x 12½"	11-1637	17.50	11-1637-12	197.50
eye					
R	3" x 9"	11-1638	15.00	11-1638-12	170.00



N



M




O




L



Q



P



R