

Hydrocollator® moist heat packs and covers

Hydrocollator® moist heat packs



- provides up to 30 minutes of soothing, moist heat
- brown tabs on most packs for rotation in clinics
- cost-effective; provides hundreds of treatments

packs only

		each	bulk (12 each)
A	standard 10" x 12"	00-1006 16.50	00-1006-12 190.00
B	neck 5" x 24"	00-1002 19.50	00-1002-12 220.00
C	half size 5" x 12"	00-1014 13.00	00-1014-12 145.00
D	oversize 15" x 24"	00-1004 40.00	00-1004-12 432.50
E	knee-shoulder 10" x 20"	00-1012 35.00	00-1012-12 385.00
F	spine small 10" x 18"	00-1008 27.50	00-1008-12 307.50
G	spine large 10" x 24"	00-1010 35.00	00-1010-12 385.00

pack and cover sets

	each	bulk (12 each)
standard pack and cover set	00-1064 37.50	00-1064-12 422.50
neck pack and cover set	00-1062 42.50	00-1062-12 472.50

Hydrocollator® heating units and cart



- stationary heating units are small enough to be kept on the counter top or cart
- water bath keeps packs at 168°F

units with packs

		L x W x H	
M	00-2102 E-1 with 4 standard	12½" x 8" x 15"	445.00
	00-2102-2 E-1 with 3 standard, 1 neck	12½" x 8" x 15"	445.00
N	00-2102-3 E-1 with 2 standard, 2 neck	12½" x 8" x 15"	445.00
	00-2802 E-2 with 2 OS, 3 std., 1 neck	14" x 12" x 16½"	912.50

carts

O	11-1390 cart with cut out (holds E-1 unit)	325.00
	11-1391 cart (holds E-2 unit)	375.00

L x W x H

P	00-4002 2 shelf MB	18¾" x 16¾" x 19½"	242.50
Q	00-4004 3 shelf MB-T	18¾" x 16¾" x 29½"	295.00

accessories

R	11-1397 hot pack tongs	15.00
	00-4228 dial thermometer	12.50

Hydrocollator® hot pack covers



all terry covers

	each	bulk (12 each)
H	standard 20" x 24"	00-1118 35.00
I	neck 9" x 24"	00-1104 27.50
J	oversize 24" x 30"	00-1126 47.50

foam-filled terry covers

	each	bulk (12 each)
K	standard 20" x 24"	00-1102 27.50
	neck 9" x 24"	00-1120 24.50
L	oversize 24" x 30"	00-1124 35.00
	standard with pocket	00-1108 22.50
	oversize with pocket	00-1110 37.50

00-4016	wall-mounted towel rack (33" x 2" with six 13" hooks)	67.50
---------	---	-------

Hydrocollator® mobile heating units



- units are large enough for a complete department
- no plumbing required; simply fill with water and plug in
- thermostatically maintains a temperature of 168°F

units with packs

		L x W x H	
00-2302	SS-2 with 8 standard	20¾" x 15½" x 33"	1,450.00
00-2402	M-2 with 12 standard	26½" x 15½" x 33½"	1,635.00
00-2502	M-4 with 24 standard	20¾" x 20" x 33"	2,975.00

accessories available with other heat pack combinations

00-4010	optional side rack	510.00
00-4230	extra shelf for side rack	47.50