


hot therapy products

TheraTemp® moist heat therapy



- place packs in microwave to heat; no water needed!
- produces moist heat by drawing moisture from the air; retains heat in the therapeutic beads
- elastic belts with hook/loop closure hold wrap in place

universal wrap packs with elastic 3 x 27" belt

A	11-1282	6" x 12"	27.50
B	11-1281	9" x 12"	35.00
C	11-1283	12" x 16"	40.00
D	11-1288	14" x 18"	50.00

specialized contour wrap packs with elastic belt

E	11-1280	cervical, 6" x 24" w/ 3 x 27" belt, 2" x 8" strip	40.00
F	11-1285	lumbar, 9" x 18" w/ 3" x 27" belt	40.00
G	11-1286	eye/sinus, 10" x 3.5" w/ 3" x 14" strap	27.50
H	11-1287	wrist, 5.5" x 12"	30.00

economy heat pads



- 3 heat settings
- washable cover
- add water to the sponge for moist heat; remove sponge for dry heat

dry heat pad

11-1130	small	12" x 15"	27.50
11-1131	large	12" x 24"	37.50

moist or dry heat pad

11-1132	small	12" x 15"	30.00
11-1133	large	12" x 24"	37.50

moist heat pads - analog and digital



- electric pad provides moist heat without adding water (fomentation)
- analog hand controller adjusts temperature
- digital hand controller adjusts temperature and treatment time
- includes flannel cover

analog pads

11-1120	king	13" x 27"	60.00
11-1121	medium	13" x 13" (pictured)	50.00
11-1123	small	4" x 14"	45.00

digital pads

11-1135	king	26" x 14"	72.50
11-1136	medium	18" x 14"	60.00
11-1137	mini	19" x 7"	50.00
11-1138	shoulder	16" x 12" (pictured)	72.50

Relief Pak® Hot Button® compress





instant treatment



- instant heat therapy without hot water or electricity
- great to use on-the-go to warm up muscles and relieve pain



- activate compress by gently flexing disk within the compress
- instantly heats to approximately 130°F (54°C)
- generates heat for at least 30 minutes
- recharge compress by boiling pad
- pad can be recharged at least 100 times



each

bulk (12 each)

I	small	3½" x 5½"	11-1025	7.50	11-1025-12	85.00
J	standard	5" x 15"	11-1026	12.50	11-1026-12	142.50
K	neck	11" x 17"	11-1027	17.50	11-1027-12	200.00
L	oversize	10" x 11"	11-1028	15.00	11-1028-12	170.00

hot water bottle



- hospital-grade rubber hot water bottle with screw stopper
- 2-quart capacity

11-1140	each	7.50
11-1140-12	dozen	82.50