

Chattanooga® Fluidotherapy® dry heat

- increase circulation and ROM with the soothing dry heat of the Fluidotherapy® unit...without having to deal with the messy cleanup and wetness associated with other therapies
- programmable; standard with automatic timer, wake-up preheater and locking casters
- adjustable controls include: temperature (maximum 125°), airflow (velocity of particulate medium), pulse mode (pressure bursts)



single FT110

- single extremity unit - for hand, wrist, elbow, ankle or foot
- includes 20 lb of Cellex™ medium
- 34"L x 11½"W x 33"H



double ULTRA 115

- treats two extremities
- removable from stand for easier lower extremity applications
- includes 40 lb of Cellex™ medium
- 34"L x 18½"W x 33"H



cellex medium

- can be used with both Fluidotherapy® and ThermoTherapy® systems
- 10 lb

FLU110D single FT110 5,850.00

FLU115D double ULTRA 115 7,400.00

MED0001 Cellex™ - 10 lbs 75.00

Thermo-therapy® dry heat

- granulated medium is kept in suspension by heated air
- thermostatically controlled heater can quickly heat medium to 125°F; temperature indicated on the unit
- built-in electric timer to control treatment time
- variable stream blowers (90-1720 RPM) control heated air streams
- standard with granulated medium and instructions

single hand TT101

- 10-12 lb capacity
- includes 10 lb medium
- 15"H x 19"D x 11"W

11-1900 single TT101 3,130.00

single hand TT101L

- 12 lb capacity
- includes 10 lb medium
- 15"H x 22"D x 11"W

11-1901 single TT101L 3,755.00

dual hand TT202

- 40-50 lb capacity
- includes 40 lb medium
- 39"H x 29"D x 20"W

11-1911 double TT202 5,270.00

stand for TT101 and TT101L

11-1920 stand for TT101 372.50
11-1921 stand for TT101L 245.00

