

electrodes

Axelgaard ValuTrove® X


- features an oversized border for easy removal
- Multi-Stick® gel adhesive provides long life and extended electrode wear
- 4 electrodes/sheet; 10 sheets/case



13-1256-10	2" round, cloth	40/cs	47.50
13-1257-10	2" square, cloth	40/cs	47.50
13-1258-10	2" round, foam	40/cs	47.50
13-1259-10	2" square, foam	40/cs	47.50

EZ-Trode electrode

- long-lasting, self-adhesive electrodes in a re-sealable pouch to preserve freshness
- thin silver layer provides superior conductivity
- soft rubber connectors
- 4 electrodes/sheet; 10 sheets/case




13-3080-10	2" round	40/cs	95.00
13-3081-10	2 3/4" round	40/cs	135.00
13-3082-10	2" x 5" rectangle*	20/cs	95.00
13-3083-10	2" square	40/cs	95.00

* 2/sheets; 10 sheets/case

V-Trode electrode


- soft rubber lead wires make these self-adhesive electrodes easy to attach to device lead wires
- long-lasting gel with a thin carbon film layer provides strong adhesion to the skin
- 4 electrodes/sheet; 10 sheets/case



13-3220-10	2" round	40/cs	60.00
13-3221-10	2 3/4" round	40/cs	82.50
13-3222-10	2" x 4" oval	40/cs	75.00
13-3223-10	2" square	40/cs	60.00

AdvanTrove® Essential electrode

- dual-layer design of adhesive hydrogel enhances performance and patient comfort and maximizes reusability
- provides uniform current throughout entire electrode for greatest treatment surface area
- 4 electrodes/sheet; 10 sheets/case



tan tricot

13-1250-10	2" square	40/cs	45.00
13-1251-10	2" round	40/cs	45.00
13-1254-10	1 3/8" round	40/cs	42.50
13-1252-10	1 3/4" x 3 3/4" rectangle	40/cs	67.50
13-1253-10	2" x 4" oval	40/cs	77.50

white foam

13-1260-10	2" square	40/cs	45.00
13-1261-10	2" round	40/cs	45.00
13-1262-10	1 3/4" x 3 3/4" rectangle	40/cs	67.50
13-1263-10	3" x 5" oval	40/cs	72.50

AdvanTrove® Essential electrode

- general use, cost-effective electrodes that provide comfort and reusability
- uses industry leading hydrogel for durability
- copper-stranded wire for enhanced current dispersion; tab and molded grip prevent wire pull out; spunlace topcoat for comfort and flexibility
- 4 electrodes/sheet; 10 sheets/case
- no expiration date



13-1240-10	2" square	40/cs	40.00
13-1241-10	2" round	40/cs	40.00
13-1242-10	1 3/4" x 3 3/4" rectangle	40/cs	55.00

Freedom-Stim™ dysphagia electrodes

- for placement on the anterior portion of the neck
- silver-coated carbon film for reduced impedance, greater conductivity and uniform current dispersion
- high quality and cost-effective (compare to Vital Stim)
- no expiration date
- dual layer of gel will not separate from carbon and is formulated for maximum reusability

A 8" round with lead, tan	13-1270-12	12 ea	170.00	13-1270-50	50 ea	700.00	13-1270-100	100 ea	1,330.00
B 8" round with snap, tan	13-1271-12	12 ea	170.00	13-1271-50	50 ea	700.00	13-1271-100	100 ea	1,330.00
C 2" x 3 1/2" butterfly with lead, tan	13-1272-12	12 ea	192.50	13-1272-50	50 ea	800.00	13-1272-100	100 ea	1,520.00
D 2" x 3 1/2" butterfly with snap, tan	13-1273-12	12 ea	192.50	13-1273-50	50 ea	800.00	13-1273-100	100 ea	1,520.00



Fabrication Enterprises, PO Box 1500, White Plains, New York 10602 USA

tel: 800-431-2830 • 914-345-9300 Orders@FabEnt.com fax: 800-634-5370 • 914-345-9800

2015-2016