

portable TENS, IF, NMES, Galvanic and EMS

Intellect® digital TENS

- battery-operated, portable dual channel digital TENS stimulation unit
- adjustable frequency and pulse duration and timer; 5 modes: burst, normal, modulated rate, width strength duration and B-Pulse; easy-to-read LCD display with convenient digital controls



07-7712 digital TENS 90.00

Intellect® economy TENS

- dual channel modified square waveform with adjustable amplitude (0-80mA), pulse frequency (2-150 Hz) and pulse width (60-25 µS) with timer
- includes 9V battery, Dura-Stick II self-adhesive electrodes and carrying case



07-7710 TENS economy 82.50

Intellect® standard TENS

- dual channel with adjustable parameters including frequency, pulse duration and 3 modulation modes and timer
- includes 9V battery, Dura-Stick II self-adhesive electrodes and carrying case



07-7600 TENS standard 90.00

accessories for Intellect® economy and standard TENS

07-7610 battery recharger 18.50
07-7609 rechargeable NiCd 9V battery 17.50

Intellect® IFC interferential

- battery-operated dual channel IFC stim
- for 2 or 4 electrode use; adjustable pulse amplitude, interference frequency and difference frequency parameters; 4 stim modes: continuous, 1/1, 8/8, and 10/10 ramped; 120V AC adapter



07-7719 IFC 240.00

Intellect® NMES

- dual channel, battery operated portable NMES stimulation unit
- includes adjustable pulse amplitude, pulse frequency and contraction/ relation time parameters



07-7715 NMES 100.00

Intellect® digital NMES


- dual channel, battery operated digital NMES stim
- includes adjustable pulse amplitude, width, frequency and contraction/relaxation time; 3 stimulation modes: synchronous, constant and alternate; easy read LCD display and convenient digital controls



07-7717 digital NMES 145.00

Mettler® TENS Sys*Stim

- 2 independently controlled channels deliver burst or modulated Symmetric Bi-Phasic square pulse wave output



A 13-1317 Sys*Stim 210 40.00
B 13-1318 Sys*Stim 210T 45.00 with timer
C 13-1314 Sys*Stim 212 85.00 with digital display

portable galvanic stimulator


- one-handed operation: current adjustment, active probe positioning and on-off control
- EMS-1C produces "Interrupted DC" (galvanic) stimulation
- EMS-2C produces "Interrupted DC and Pulsed" (faradic) stimulation
- includes a 9V battery



13-1310 EMS 2 297.50
13-1311 EMS 1 250.00

mini wireless TENS/EMS

- single channel output, 3V CR2025 battery (included), output 0~90 mA
- timer: 30 min. backward to auto shut off
- dual phase square wave
- LCD display with timer/mode/density
- modes: 5 TENS programs; 5 EMS fitness programs; contraction (on time) 2-4 sec; relaxation (off time) 2-4 sec



13-1360 mini unit 25.00

Chattanooga® Primera TENS/NMES

- extensive range of programs – 8 TENS and 6 NMES programs
- TENS/NMES with HAN Waveform. Low (2Hz) and high (15 or 70Hz) frequencies are sequenced to better manage pain
- complete package includes: unit, carrying case, battery, lead wires, Dura-Stick electrodes, back support belt and user's manual



07-7615 Primera 147.50