

## FabStim<sup>®</sup> TENS 2-channel

- 3 modes (Burst, Modulation, Constant) and timer
- adjustable pulse rate and width
- includes carry case, stimulator, lead wires, 4 self-stick electrodes and 9V battery



13-1323 analog TENS 32.50  
13-1313 Dual Channel TENS 37.50

## digital EMS/TENS 4-channel


- 5 TENS modes, 3 EMS modes
- variable pulse rate and width
- variable ramp time
- variable contraction/relaxation time
- waveform: Asymmetrical Bi-Phasic square pulse
- 2 timers and compliance meter
- comes with 8 electrodes



13-1349 digital EMS/TENS 147.50

## IF 2-channel w/ AC adapter

- 4 modes: continuous, 1/1 abrupt, 8/8 abrupt, 10/10 ramped
- frequency shift: from 25% below set frequency to 55% above
- carrier frequency: 4000Hz
- modulating frequency: 1-150Hz
- wave: Symmetrical Bi-Phasic square pulse
- includes stimulator, lead wires, 4 self stick disposable electrodes, 9V battery and AC adapter



13-1303 interferential 75.00

## FabStim<sup>®</sup> digital TENS 2-channel

- 5 modes
- adjustable pulse rate and width
- waveform: Asymmetrical Bi-Phasic square pulse
- timer (5-60 mins or continuous) and compliance meter
- includes carry case, stimulator, lead wires, 4 self stick electrodes and 9V battery



13-1308 digital TENS 40.00

## microcurrent TENS 2-channel

- dual channel microcurrent / micro TENS
- 3 Frequencies: 0.3Hz, 8Hz, 80Hz
- changes polarity each second
- waveform: Modified square Bi-Phasic
- timer: constant, 20 min., 40 min.
- includes carry case, stimulator, lead wires, 4 self stick electrodes, 9V battery and instructions



13-1322 microcurrent TENS 60.00

## EMS 2-channel

- variable pulse: rate 5, 30, 100(Hz), width 250(μs)
- variable ramp time: 1,3 or 5 seconds
- variable contraction/relaxation time
- wave: 3-types
- includes carry case, stimulator, lead wires, 4 self stick electrodes, 9V battery and instructions



13-1306 EMS 40.00

## EMS/TENS 2-channel

- 3 TENS modes, 2 EMS modes
- variable pulse rate and width
- variable ramp time
- variable contraction/relaxation time
- waveform: Asymmetrical Bi-Phasic square pulse
- timer and compliance meter
- includes carry case, stimulator, lead wires, 4 self stick electrodes and 9V battery



13-1348 EMS/TENS 60.00

## IF 4K digital interferential

- digital interferential with multiple treatment modes
- AC adaptor provides up to 2 hours of treatment



13-1379 IF 4K interferential 130.00

## EMS 2-channel digital

- 3 modes
- adjustable pulse rate and width
- variable ramp time: 1, 3 or 5 seconds
- variable contraction/relaxation time; timer and compliance meter
- wave: 3-types
- includes carry case, stimulator, lead wires, 4 self stick electrodes, 9V battery and instructions



13-1309 digital EMS 37.50