

Chattanooga® Intellect® LEGEND XT™ Therapy Systems

Intellect® LEGEND XT™ Therapy Systems

Electronic "platform" is capable of providing multiple modalities in one therapy station. You can start simple (stim and ultrasound) and add modalities (additional stim, laser, battery module and patient management software).


- stim unit offers 6 clinical waveforms including: interferential, Russian, bi-phasic, high volt, microcurrent and premodulated
- ultrasound delivers 1 and 3.3 MHz through a 5cm² applicator. Multiple waveforms available in combination mode
- built-in patient data card device can record patient modality and pain profile information
- user-defined protocols customize any waveform to your specific need
- optional mobile cart has 6 storage bins to minimize cord and cable hassle

- each unit ships with 5 Patient Data Cards, 4 each 3" inch diameter carbon electrodes (2 red, 2 black), 4 each Dura-Stick® II 2 3/4" diameter electrodes and 1 2 1/2" x 24" Nylatex® wrap

MAP (minimum advertised price)



stim unit (left)
combo unit: stim plus
ultrasound (right)
both shown with
optional mobile cart



stim units

00-2763	2-channel stim	MAP	2,870.00
00-2793	2-channel stim with cart	MAP	3,350.00
00-2786	4-channel stim	MAP	3,620.00
00-2797	4-channel stim with cart	MAP	4,075.00
00-2780	mobile cart		457.50
02-7508	interrupt switch channel 1,2		115.00
02-7079	interrupt switch channel 3,4		115.00

combo units (stim plus ultrasound)

standard with 5cm² (1 and 3.33 MHz) applicator, 9 oz conductor gel

00-2760	2-channel combo	MAP	4,100.00
00-2791	2-channel combo with cart	MAP	4,585.00
00-2788	4-channel combo	MAP	4,800.00
00-2795	4-channel combo with cart	MAP	5,260.00
07-2852	rubber 3" carbon electrode black		22.00
07-2853	rubber 3" carbon electrode red		22.00