

Chattanooga® Intellect® TRANSPORT™ therapy systems

Transport™ 2-channel stimulation

- 2-channel output with 4 standard waveforms: interferential, premodulated, HV and Russian
- logical and simple user interface with 10 programmable memory positions
- built-in positional stand allows desktop or wall-mount use
- standard with four ¾" diameter Dura-Stick® II electrodes, four 3" diameter carbon electrodes (2 red, 2 black), lead wires and one 2½" x 24" Nylatex® wrap

00-2783K 2-channel stimulation unit
 00-2894K stim with mobile cart
 00-2783B stim with carry bag and battery pack

MAP 2,415.00
 MAP 2,965.00
 MAP 2,890.00

MAP
 (minimum advertised
 price applies)



Transport™
stim unit

Transport™ 2-channel stimulation / ultrasound combination

- electrotherapy and ultrasound work seamlessly to heat and stimulate muscles for more rapid healing and better patient outcome
- features Transport™ 2-channel stimulation and ultrasound

00-2738K combo: 2-channel stim and ultrasound w/ 5cm² applicator
 00-2738KC combo with mobile cart
 00-2738KB combo with bag and battery pack

MAP 3,350.00
 MAP 3,905.00
 MAP 3,830.00



Transport™
combo unit

Transport™ ultrasound

Transport™
ultrasound unit



- Transport™ ultrasound delivers 1 and 3.3 MHz frequency in continuous (100%) or pulsed (10%, 20%, 50%) format
- includes 10 programmable memory positions
- output is a clean, culminated beam with a BNR that is of the lowest available (max 5:1)
- built-in positional stand facilitates desktop or wall-mount use
- powered by battery or 110V AC
- standard with one 9 oz ultrasound gel


00-2782K ultrasound w/ 5cm² applicator
 00-2893K ultrasound with cart

MAP 1,770.00
 MAP 2,325.00

mobile cart

- optional mobile cart and adapter for Intellect® Transport™ unit

00-2780 Legend XT™ 457.50
 00-2891K Transport™ cart+adapter 560.00
 00-2884 adapter only 105.00



cart
adapter



cart with
adapter,
stim/ultrasound
sold separately

Transport™ battery pack



- nickel metal hydride battery stores in bottom of Intellect® Transport™ units

02-7478 battery 340.00
 00-2911K battery & bag 480.00

Transport™ carry bag



- ergonomically designed bag for easy transport of Transport™ models
- holds unit securely in place and fits accessories

02-7467 carrying bag 142.50

additional soundheads



- additional applicators for Legend XT™ (on page 112) and Transport™ ultrasound

02-7381 1cm² 360.00
 02-7382 2cm² 360.00
 02-7383 5cm² 360.00
 02-7384 10cm² 360.00