

traction devices

HipTrac® hip traction

- use for hip osteoarthritis, FAI/labral pathology, joint stiffness, pain or muscle tightness
- use in patients' home without assistance of a second person
- use in a clinic to prepare patient for joint and soft mobilization
- reduce pain and increase mobility to improve exercise compliance and performance
- adjustable ankle binding

50-0920 HipTrac™ 895.00



Posture-Pro™ cervical traction

pro 1000 cervical
disc hydrator



pro 2000 lumbar
elliptical back rocker



pro 2500 Penta Vec



- up/down motion promotes disc lubrication/joint nutrition
- portable, easy to operate

50-1040	Posture Pro™ 1000 cervical disc hydrator	200.00
50-1041	Posture Pro™ 2000 lumbar elliptical back rocker	200.00
50-1042	Posture Pro™ 2500 Penta Vec	280.00
50-1045	Posture Pro™ 4100 full spine combo	375.00

pro 4100 full spine combo



Pronex® pneumatic cervical traction

- even distraction in the anterior and posterior cervical discs
- patient controls the amount of traction applied
- squeezing an inflator bulb increases the pressure, release knob gently reduces the pressure

50-1050	regular (14"-16" neck)	395.00
50-1051	large (16"-18" neck)	395.00
50-1052	wide (18"+ neck)	395.00

Pronex® cervical traction

