

Chattanooga® continuous passive motion (CPM)

OptiFlex-K1™ knee

2
year
warranty



features:

- weighs just 23 lb
- only 38" (97 cm) long
- wide low-profile, closed-base frame for stability
- faster motor than the K3
- accommodates patients 4.0 - 6.6 ft (1.2 - 2 m)
- reversible footplate for short leg length
- on-board power supply, ROM: -10° to 120°

knee CPM unit

00-2025	CPM with classic hand control	2,805.00
00-2026	CPM with standard hand control	2,805.00
00-2028	CPM with comfort hand control	2,805.00

hand control only

03-7413	classic hand control only	527.50
03-7411	standard hand control only	527.50
03-7412	comfort hand control only	527.50

patient kits

02-0533	OptiFlex® patient kit	27.50
01-1432	Breg patient kit	27.50

straps

04-0128	OptiFlex® 34" strap kit	8.50
---------	-------------------------	------

OptiFlex-K1 offers 3 unique hand controllers capable of matching the K1's performance settings precisely to the specific needs of the patient. This assures that the patient will receive the right therapy treatment.



Classic



Standard



Comfort

	Classic	Standard	Comfort
Key lock function	✓	✓	✓
Pauses	✓	✓	✓
Speed	✓	✓	✓
Timer		✓	✓
Reverse on load		✓	✓
Warm up protocol		✓	✓
Default setting		✓	✓
Transport setting		✓	✓
Chip card selectable		✓	✓
Special Functions			
EROM repeat			✓
Stretch protocols			✓
Workout protocol			✓
Comfort protocol			✓
Therapy documentation			✓

OptiFlex® knee CPM trolley

- hi-lo knee CPM trolley with multifunctioned trolley design
- easy storage and transportation
- allows for proper CPM positioning
- adjustable height from 18¼" to 28"

positioning trolley

00-4520	hi-lo knee trolley	697.50
---------	--------------------	--------

pictured with
OptiFlex3®

