

# Kinetec® continuous passive motion (CPM)

## Breva ankle



- provides anatomically correct motion in two planes
- offers plantar flexion, dorsiflexion, and inversion and eversion patterns
- quick and easy set up and change-overs
- designed for use in a chair or bed

13-1072 Breva ankle CPM 5,575.00

## Maestra® portable hand



- short set-up time by incorporating a low profile rigid splint
- ROM: 15° hyperextension to 270° flexion (composite fist)
- lightweight and easy to operate
- anatomically shaped drive bar enhances abduction and adduction
- up to 15 minutes of pause at extension and flexion allows for a controlled stretch and rest period
- convenient display of angle

13-1165 Maestra® hand CPM 5,425.00

## Maestra® hand/wrist



- ROM:
  - fingers:
    - 45° hyperextension
    - 27° flexion (full list)
  - thumb:
    - full opposition
  - wrist:
    - -50° extension,
    - -90° flexion
    - -90° pronation
    - -90° supination

13-1168 Maestra® hand/wrist CPM 6,800.00

## Centura™ and 6080 elbow

- provides anatomically correct elbow extension/flexion in an abducted position
- quick and easy change-over from an existing shoulder pattern
- ROM: extension 10°, flexion 135°
- extension with fixed, adjustable pronation/supination in an abducted, rotated plane



6080 elbow CPM

### Centura™

13-1155 Centura™ elbow (CEM) 3,475.00  
13-1158 Centura™ patient kit 137.50

### 6060 elbow CPM

13-1074 6080 elbow 6,450.00  
13-1078 6080 patient pad 47.50

## Centura™ shoulder

- helps prevent joint stiffness, soft tissue contractures, atrophy
- patient set-up is anatomical, comfortable and quick with a color-coded set-up system
- easy change for left or right
- breaks down for convenient transportation
- includes carrying case (for detachable parts), pad kit and a reusable chair cover



13-1157 Centura™ shoulder CPM 6,725.00  
13-1151 replacement patient pad kit 40.00  
13-1153 disposable chair cover 30.00

## shoulder bed/wheelchair

- helps prevent joint stiffness, soft tissue contractures, atrophy
- patient set-up is anatomical, comfortable and quick with a color-coded set-up system
- easy change for left or right
- use while the patient lies in bed or sits in a wheelchair
- programming through the hand control allows for extension/flexion movements
- breaks down for convenient transportation
- ROM: arm in full extension: 0°-110°
- elbow flexed 90°: 30°-110°



13-1156 Centura™ bed/wheelchair 6,250.00  
13-1161 patient kit 40.00