

Kinetec® knee continuous passive motion (CPM)

R-Flex



- lightweight and portable, only 27.5 lbs with self-contained handles
- compact design, 29.6" long
- ROM: -10° hyperextension to 120° of flexion allows for treatment of all knee pathologies
- widest carriage available at 13"
- easy to read LED display
- easy calibration procedure
- foam patient pad kit included
- 14"W x 20"H x 29.6"L

13-1230	R-Flex knee machine	3,100.00
13-1231	R-Flex management system	560.00

Kinetec® Prima™ Advance / XL Bariatric

- easy to use and adjust
- simple straightforward movement adjustments for ROM limits, speed and patient sizing
- 37" length simplifies use in a hospital bed or in small spaces anywhere
- ROM: -5° to 115° of flexion
- 24 lb



standard


13-1070	Prima™ Advance	3,350.00
---------	----------------	----------

bariatric



- bariatric unit has 12½" extra-wide carriage and a 500 lb capacity
- | | | |
|---------|-----------|----------|
| 13-1071 | Prima™ XL | 4,875.00 |
|---------|-----------|----------|


Kinetec® Spectra™ / Essential



- anatomically corrects alignment to increase patient compliance and decrease joint stress
- wide carriage (10") accommodates all patients
- ROM: hyperextension (-10°) through full extension (120°)
- load reverse for patient protection
- patient lockout prevents program changes
- pause settings for 15 minutes
- 24 lb

13-1185	Spectra™ with basic hand control	4,050.00
13-1189	Spectra™ Essential with advanced hand control	4,050.00
13-1186	pad set	22.50

Kinetec® Performa™



- leg carriage for automatic motion
- adult and pediatric capability
- programmable muscle stim capability
- includes a carrying case for easy transportation
- built in therapy timer
- attached sliding hip bar
- manual set-up mode
- pad system and straps

13-1170	Performa™ knee CPM	5,025.00
13-1171	replacement patient pad set	22.50

Kinetec® Pedala Active Trainer



- the ideal complement to ensure full knee recovery
- isometric, isotonic and physiological work; muscular reinforcement
- exceptional patient compliance and autonomy
- patient self-controls the movement, the speed and the ROM
- 27 lb; 35"L x 14"W x 15"H




with seat adapter

13-1190	Kinetec® Pedala	2,950.00
13-1191	Pedala seat adapter	375.00
13-1192	exchange complete pads set	50.00

Kinetec® accessories

adjustable height trolley

- provides transport positioning and storage
- 3" caster wheels; two locking
- use with Spectra™ Performa™, Prima™ XL and R-Flex knee CPMs



13-1187	CPM knee trolley	742.50
---------	------------------	--------

Kinetec® CPM cart

- transport and store knee CPM units and accessories
- resistant to commonly used disinfectants
- 27½"L X 14½"W x 29½"H



13-1148	CPM cart	365.00
---------	----------	--------