

AP / LAL pressure pump and pad

AP / LAL alternating pressure / low air loss

Med-Aire 8



- mattress can operate in alternating or static pressure mode
- 8" deep air bladders are easily removed and replaced
- CPR valve allows for rapid deflation
- 80"L x 36"W x 8"H

system

13-0850 Med-Aire 8 770.00

accessories

13-0751 repl. pump 250.00
13-0753 repl. cover 90.00
13-0754 repl. bladder without air flow 22.50
13-0755 repl. bladder with air flow 20.00

Med-Aire Plus



- digital pump allows for customization
- 8" deep air bladders are easily replaced
- 3 diagnostic alarms
- CPR valve allows for rapid deflation
- 80"L x 36"W x 8"H

system

13-0860 Med-Aire Plus 1,925.00

accessories

13-0761 repl. pump 835.00
13-0763 repl. cover 230.00
13-0764 repl. bladder without air flow 82.50
13-0765 repl. bladder with air flow 72.50

Balanced Aire Powered



- pressure pump redistributes pressure through 3 zones
- system is sloped downward in the foot section to reduce pressure in heel zone
- 80"L x 7"H

system width

13-0870-35 35" air/foam 1,652.50
13-0870-42 42" air/foam 2,312.50
13-0870-48 48" air/foam 2,752.50
13-0870-54 54" air/foam 3,082.50

accessories

13-0781 repl. pump 112.50

Balanced Aire Non-Powered



- 3 pressure zones independently self-adjust to patient's need
- system is sloped downward in the foot section to reduce pressure in heel zone
- 80"L x 7"H

system width

13-0880-35 35" air/foam 1,540.00
13-0880-42 42" air/foam 2,200.00
13-0880-48 48" air/foam 2,640.00
13-0880-54 54" air/foam 2,970.00

Altadyne™ 750000



- pressure control knob allows firmness of mattress to be adjusted
- LED indicator lights for low pressure alert
- lighted ON/OFF indicator switch
- CPR valve for rapid deflation

system

13-0800 Altadyne™ 1,185.00

accessories

13-0701 repl. pump 405.00
13-0702 mattress only 860.00
13-0703 mattress cover 375.00
13-0705 repl. cell 30.00
13-0706 foam base 160.00

Altadyne™ Plus



- touch pad adjustment with LED indicators
- 10, 15, 20, and 25 minute cycle times
- lockable control panel
- low pressure and power failure LED indicators and audible alarms
- CPR valve for rapid deflation

system

13-0810 Altadyne™ Plus 1,300.00

accessories

13-0711 repl. pump 815.00
13-0712 mattress only 1,605.00
13-0713 mattress cover 255.00
13-0714 cell with air loss 90.00
13-0715 cell without air loss 50.00
13-0716 foam base 160.00

Altadyne™ Sensor Plus



- alternating pressure relief support surface Auto-Adjusting technology
- 10 minute cycle time
- lockable control panel
- low pressure and power failure LED indicators and audible alarms
- CPR valve for rapid deflation

system

13-0820 Sensor Plus 1,500.00

accessories

13-0721 repl. pump 840.00
13-0722 mattress only 1,500.00

Altadyne™ XD



- suitable for use in the prevention and treatment of stages 1 - 4 pressure ulcers
- cycle time adjusts from 5-20 minutes, in 5 minute increments
- provides 6 options for mattress pressure, plus an optional static mode

system

13-0830 Altadyne™ XD 4,950.00

accessories

13-0731 repl. pump 2,025.00
13-0732 mattress only 4,385.00
13-0733 mattress cover 210.00
13-0735 cell w/o air loss 150.00