

AP pressure pump and pad

APP alternating pressure pad and pump system

Med-Aire Alternating Pressure Pump and Pad System



variable pressure pump

- 130 individual bubble cells for maximum comfort and effectiveness
- variable pressure pump, 5-min. cycle time
- variable pressure setting on deluxe pump allows comfort setting for maximum comfort and compliance
- durable, heavy gauge, premium vinyl bubble pad provides superior therapeutic treatment
- built-in bracket for easy mounting to bed

system

| | | |
|---------|---|--------|
| 13-0950 | variable pressure pump and pad | 127.50 |
| 13-0951 | variable pressure pump and pad with end flaps | 135.00 |

accessories

| | | |
|---------|--|--------|
| 13-0955 | variable pressure pump only | 112.50 |
| 13-0956 | standard pad, 78"L x 36"W x 2.5"H | 37.50 |
| 13-0957 | pad with end flaps, 123"L x 36"W x 2.5"H | 42.50 |

Med-Aire Fixed Pressure Pump and Pad System



fixed pressure pump

- 130 individual bubble cells for maximum comfort and effectiveness
- fixed pressure pump
- durable, heavy gauge, premium vinyl bubble pad provides superior therapeutic treatment
- built-in bracket for easy mounting to bed

system

| | | |
|---------|-----------------------------|--------|
| 13-0960 | fixed pressure pump and pad | 120.00 |
|---------|-----------------------------|--------|

accessories

| | | |
|---------|-----------------------------------|--------|
| 13-0965 | fixed pressure pump only | 100.00 |
| 13-0966 | standard pad, 78"L x 36"W x 2.5"H | 37.50 |

AQ1000 alternating pressure system

bubble pad



standard alternating pressure pump



- bubble style pad effectively reduces pressure
- pump mounts to bed with built-in mounting bracket; features an illuminated on/off switch
- standard alternating pressure pump with 5-min. cycle
- 78"L x 34"W x 2 1/2"H
- 275 lb maximum weight capacity

| | | |
|---------|------------------------------|--------|
| 13-0900 | standard pump and bubble pad | 165.00 |
| 13-0901 | standard pump only | 80.00 |
| 13-0991 | replacement pad | 35.00 |

AQ2000 deluxe variable pressure system

bubble pad



deluxe variable pump



- bubble style pad effectively reduces pressure
- pump mounts to bed with built-in mounting bracket; features an illuminated on/off switch
- deluxe variable pressure pump with 5-min. cycle
- 78"L x 34"W x 2 1/2"H
- 275 lb maximum weight capacity

| | | |
|---------|-----------------------------------|--------|
| 13-0910 | deluxe variable pump & bubble pad | 175.00 |
| 13-0911 | deluxe variable pump only | 105.00 |
| 13-0991 | replacement pad | 35.00 |