



- Game Ready® with active compression delivers cold therapy that works faster, penetrates deeper and lasts longer than other treatment modalities
- the cold compression wraps incorporate an active temperature exchange unit that maintains cool temperature and a pneumatic compression
- the system has a control unit and full range of ergonomic wraps for all body parts

control unit		
13-2500	GRPro 2.1 control unit includes AC adapter and 6 ft connector hose	2,830.00

wraps		additional sleeves		
upper extremity		upper extremity		
13-2519	430.00	13-2561	60.00	hand/wrist
13-2520	230.00	13-2562	20.00	hand with sleeve
13-2521	200.00	13-2563	10.00	hand with disposable coverlet*
13-2522	320.00		not available	wrist (one size)
13-2523	350.00	13-2564	30.00	elbow (one size)
13-2524	490.00	13-2565	50.00	flexed elbow*
13-2525	490.00	13-2566	50.00	shoulder right, medium (33-45" chest)
13-2526	490.00	13-2567	50.00	shoulder right, large (40-55" chest)
13-2527	490.00	13-2568	50.00	shoulder left, medium (33-45" chest)
13-2528	490.00	13-2569	50.00	shoulder left, large (40-55" chest)
lower extremity		lower extremity		
13-2510	470.00	13-2552	60.00	below knee, traumatic amputee, large
13-2511	420.00	13-2553	30.00	ankle, large (men's shoe sizes up to 11)
13-2512	420.00	13-2554	30.00	ankle, x-large (men's shoe sizes 12-18)
13-2513	490.00	13-2555	50.00	half leg boot, large
13-2514	430.00	13-2556	30.00	knee, straight (one size)
13-2515	500.00	13-2557	30.00	knee, articulated (one size)
mid body		mid body		
13-2516	490.00	13-2558	50.00	hip/groin right
13-2517	490.00	13-2559	50.00	hip/groin left
13-2518	490.00	13-2560	50.00	back

accessories			
13-2540	carry bag	190.00	13-2541 bag for 10 wraps 100.00
13-2544	rechargeable battery pack	320.00	13-2542 5 ft dual connector hose 440.00
13-2545	battery pack cable	20.00	13-2543 6 ft connector hose 200.00

*does not include heat exchanger

PresSsion® sequential / gradients sequential compression therapy

4 chamber sequential 651-4

- intermittent gradient sequential
- 180 second air eject after shut off
- inflation (15 sec/chamber), deflation (12 sec/cycle), timer (30-180 min), pressure (20-80 mmHG)



pump only

00-4307	4 chamber	750.00
---------	-----------	--------

4 chamber sequential / gradient sleeves

arm		length	
04-3207	medium	17"	340.00
04-3208	large	29"	340.00
leg		length	
04-3197	small	10"	250.00
04-3198	small, wide	20"	250.00
04-3199	small, x-wide	26"	250.00
04-3200	medium	29"	250.00
04-3201	medium, wide	32"	250.00
04-3202	medium, x-wide	35"	250.00
04-3203	large	29"	250.00
04-3204	large, wide	32"	250.00
04-3205	half leg	35"	150.00
04-3206	half leg, wide	29"	150.00

blocker

04-3221	4 chamber blocker	4.00
---------	-------------------	------



medium leg

8 chamber gradient sequential 652-8

- gradient or pressure modes
- continuous operation with intermittent loading
- uncorking and focus modes
- 60 second air eject after shut off
- inflation (25-30 sec/chamber), deflation (12 sec/cycle), timer (15-180 min), pressure (20-129 mmHG)



pump only

00-4309	8 chamber	1,140.00
---------	-----------	----------

8 chamber sequential / gradient sleeves

arm		length	
04-3217	medium	17"	340.00
04-3218	large	29"	340.00
04-3219	foot/ankle	10"	340.00
04-3220	half	20"	340.00
leg		length	
04-3209	small	26"	250.00
04-3210	small, wide	29"	250.00
04-3211	small, x-wide	32"	250.00
04-3212	medium	35"	250.00
04-3213	medium, wide	29"	250.00
04-3215	large	32"	250.00
04-3216	large, wide	35"	250.00

blocker

04-3222	8 chamber blocker	6.00
---------	-------------------	------



medium leg