

## stainless steel whirlpools

- whirlpools have a single ½ horsepower motor to power the agitator, spot spray and to empty tub
- dial varies the degree of agitation by controlling the amount of air introduced into the water
- agitator can be raised and lowered to change the directional flow; can be removed from the shaft if service is needed
- comes with thermometer
- available with a polished satin finish or an electrostatically applied color
- tank can be emptied through the gravity drain or by using whirlpool turbine assembly and drainage hose

mobile slant back lo-boy



mobile lo-boy



stationary hi-boy



mobile extremity unit



mobile hi-boy



- **mobile whirlpool** is mounted on 4 heavy-duty 5" casters
- **stationary whirlpool** is designed for permanent installation with direct plumbing connections to the water supply and drain

### lo-boy whirlpool

- long and shallow units are designed for easy access
- lower height perfect for pediatric and geriatric patients
- can provide treatment to various extremity regions along hip and back areas

		L" x W" x D" x H"	gallons	
42-1214	mobile slant back SB-100-S	76 x 40 x 21 x 24	100	13,610.00
42-1215	mobile L-75	52 x 24 x 18 x 24	75	6,800.00
42-1216	mobile L-90-M	60 x 24 x 18 x 24	90	7,225.00
42-1217	mobile L-105-M	66 x 24 x 18 x 24	105	7,650.00
42-1211	stationary L-90-S	60 x 24 x 18 x 24	90	7,585.00
42-1212	stationary L-105-S	66 x 24 x 18 x 24	105	8,030.00

### hi-boy whirlpool

- provides complete coverage to the back, hip and shoulder areas

		L" x W" x D" x H"	gallons	
42-1225	mobile H-60-M	36 x 20 x 28 x 34	60	5,685.00
42-1228	mobile H-105-M	48 x 24 x 28 x 34	105	6,410.00
42-1223	stationary H-105-S	48 x 24 x 28 x 34	105	6,730.00

### extremity whirlpool

- compact; mobile units(\*) compatible with Hydrolift™ base (right)

		L" x W" x D" x H"	gallons	
42-1245*	mobile podiatry P-10-M	22 x 13 x 12 x 18	10 *	4,085.00
42-1255*	mobile foot/ankle/arm E-15-MU	25 x 13 x 15 x 28	15 *	5,060.00
42-1256*	mobile foot/ankle/calf/arm E-22-MU	28 x 15 x 18 x 34	22 *	5,215.00
42-1258	mobile foot/ankle/calf/knee/arm E-45-M	32 x 15 x 25 x 31	45	4,890.00
42-1251	stationary foot/ankle/calf/arm E-22-SP	28 x 15 x 18 x 34	22	5,400.00
42-1253	stationary foot/ankle/calf/knee/arm E-45-S	32 x 15 x 25 x 31	45	5,525.00



Hydrolift™ base with mobile extremity unit

### Hydrolift™ whirlpool base

- Hydrolift™ base, when attached to a whirlpool (as above), is essentially two whirlpools in one allowing you to treat lower and upper extremities with the same unit
- raises the mobile extremity and podiatry whirlpool up to 3 feet off the floor
- turn-handle raises the whirlpool quickly and easily, even when completely filled

#### Hydrolift™ base alone

42-1060	for 10-gal unit	2,485.00
42-1061	for 15-gal unit	2,755.00
42-1062	for 22-gal unit	2,865.00

#### Hydrolift™ base with mobile whirlpool

42-1070	with 10-gal unit	6,290.00
42-1071	with 15-gal unit	6,945.00
42-1072	with 22-gal unit	7,545.00