


height adjustable folding walkers



- 1" diameter aluminum construction
- U-shape frame creates greater clearance
- folds for easy transport or storage
- available with rubber tips or front wheels and rear glider tips
- height adjustable: adult 32-39"; junior 25-32"; oversize 32-39"

walker	weight capacity (lb)		each	case - qty
adult, no wheels	350	43-2100	57.50	43-2100-4 210.00
junior, no wheels	350	43-2102	60.00	43-2102-4 215.00
oversize, no wheels	500	43-2104	110.00	43-2104-2 215.00
adult, 3" wheels	350	43-2110	70.00	43-2110-4 247.50
adult, 5" wheels	350	43-2120	65.00	43-2120-4 230.00
junior, 5" wheels	350	43-2122	70.00	43-2122-4 247.50

walker platform attachment



- provides the ability to bear weight on the arm and relieve pressure for those with hand injuries or arthritis
- lightweight aluminum
- no tools needed to attach to the right or left side of walker
- may be mounted on inside or outside of walker frame

43-2130 platform attachment 115.00

Dycem® strips

- adhesive backed Dycem® strip secures onto any surface for non-slip grip
- 16"L x 1½"W strips, 3 each
- specify: blue (B), yellow (Y)



walker balls



- protects floor from scuff marks
- pre-cut to fit most walkers

50-1560 B/Y strips (3) 10.00 43-2136 balls (2 each) 20.00

rollators



shown with basket

4-wheel



3-wheel

- 8" casters for indoor/outdoor use
- folds for easy transport or storage
- removable, foldable backrest
- 3-wheel comes standard with basket, tray and pouch
- 4-wheel model has a padded seat for resting and is available with push or loop braking system

rollator	weight capacity (lb)	with loop brake	with push brake
4-wheel, blue	300	43-2150 185.00	43-2140 210.00
4-wheel, red	300	43-2152 185.00	
3-wheel, blue	300	43-2160 160.00	
3-wheel, red	300	43-2164 160.00	
optional basket			
4-wheel rollator	43-2170	32.50	
3-wheel rollator	43-2171	40.00	

walker rest seat



- comfortable, fold-down padded seat
- 2" foam-padded vinyl cover cleans easily
- two hook and loop fastener straps secure the seat to the upright or horizontal rails

43-2180 rest seat 37.50

adjustable glides



- for use with all 1" dia folding walkers
- durable composite construction ensures long life

43-2135 glides (2 each) 17.00

reacher clip

- keeps reachers attached to walker or wheelchair frame
- larger side of double clip snaps onto wheelchair, walker frame or bed rail; smaller side snaps onto the reacher



43-2181 reacher clip 10.00