

wheelchairs and accessories

transport chair



- 8" composite wheels
- rear wheel locks
- light, strong and durable aluminum frame
- swing-away footrests
- back folds for storage and transport
- seatbelt for added safety
- carry pocket on backrest
- blue frame, black upholstery
- 300 lb capacity; chair weight: 19 lb

| | width | |
|---------|----------|--------|
| 43-2200 | 17" seat | 275.00 |
| 43-2210 | 19" seat | 330.00 |

FabLife™ transfer boards



- tapered boards for those with lower extremity disabilities
- 350 lb capacity

solid wood

| | | |
|---------|--------------------|-------|
| 50-3000 | 8 x 24" no cutouts | 30.00 |
| 50-3001 | 8 x 30" no cutouts | 35.00 |
| 50-3004 | 8 x 24" 2 cutouts | 35.00 |
| 50-3005 | 8 x 30" 2 cutouts | 40.00 |

plastic

| | | |
|---------|--------------------|-------|
| 50-3007 | 8 x 28" no cutouts | 45.00 |
|---------|--------------------|-------|

dual axle wheelchair



- powder-coated steel frame is lightweight, strong and durable
- padded armrests, plastic footplates
- vinyl upholstery is durable and easy to clean
- urethane tires on composite wheels; push-to-lock brakes
- 8" front casters, adjustable in 3 positions
- carry pocket on backrest provides additional convenience
- swing-away footrests or swing-away elevating legrests with adjustable length riggings
- fixed or detachable desk arms
- 250 lb capacity

fixed arm

| swing-away footrest | | width | swing-away elevating legrest | |
|---------------------|--------|----------|------------------------------|--------|
| 43-2221 | 290.00 | 16" seat | 43-2222 | 317.50 |
| 43-2251 | 290.00 | 18" seat | 43-2252 | 317.50 |

detachable arm

| swing-away footrest | | width | swing-away elevating legrest | |
|---------------------|--------|----------|------------------------------|--------|
| 43-2230 | 317.50 | 16" seat | 43-2231 | 340.00 |
| 43-2260 | 317.50 | 18" seat | 43-2261 | 340.00 |

oxygen tank holders

- attach to handles and frame; tank position does not interfere with user's head or back
- heavy duty, waterproof nylon

mini tank holder



- 2 pockets hold 2 M6 tanks or supplies

43-2280 mini holder 25.00

standard tank holder



- fits size "D" and "E" oxygen tanks

43-2281 standard holder 22.50


Skillbuilders® arm elevator



- fits into all standard wheelchair arm support trays
- coated arm elevator can be kept sanitary

30-1410 arm elevator 65.00

Skillbuilders® knee spreader



- comfortably maintains 4" of knee separation to prevent hip abduction
- hook and loop straps hold wedge in place
- coated knee spreader can be kept sanitary
- impervious to liquid

30-1411 knee spreader 45.00