

wheelchair cushions and accessories

FabLife™ wheelchair trays

- universal design fits all standard wheelchairs
- raised rim, cut-out design; attaches with Velcro® straps



sand color plastic tray



wood finish

50-1300	economy hardboard	15.00
	24"W x 20"D x 1/2"H	
50-1301	standard walnut laminate, 24"W x 20"D x 1/2"H	35.00

clear acrylic

50-1302	with rim, straps	55.00
	23"W x 19"D x 1/4"H	

plastic

43-2293	gray plastic for 16"-20" seat widths, 24"W x 20"D x 1/2"H	65.00
43-2294	sand color, contoured, round edges, 2 cup insets, 24"W x 21"D	170.00

wheelchair desk



- ideal for laptop computer users
- storage area below lid
- fits 16"-20" standard arm wheelchairs
- 24"W x 22"D x 3 1/2"H

43-2290	laptop desk	100.00
---------	-------------	--------

clamp-on drink holder



- slotted on both sides to hold a mug or cup with handles
- fits bottles and cans
- clamp has easy-to-tighten lever

43-2286	drink holder	10.00
---------	--------------	-------

Amara® wheelchair cushions with removable anti-microbial black cover

Amara 100



- for low/moderate risk of skin breakdown
- minimal postural support

W x D x H

50-1800	16 x 16 x 3"	67.50
50-1801	18 x 18 x 3"	82.50
50-1802	20 x 18 x 3"	85.00


Amara 400

- for moderate/high risk of skin breakdown; with solid gel insert
- significant postural support

W x D x H

50-1830	16 x 16 x 3"	172.50
50-1831	18 x 18 x 3"	212.50
50-1832	20 x 18 x 3"	247.50

Amara 300



- for moderate/high risk of skin breakdown; with foam insert
- minimal postural support

W x D x H

50-1820	16 x 16 x 3"	112.50
50-1821	18 x 18 x 3"	130.00
50-1822	20 x 18 x 3"	142.50

Amara 450




- for moderate/high risk of skin breakdown; with viscous gel insert
- significant postural support

W x D x H

50-1840	16 x 16 x 3"	310.00
50-1841	18 x 18 x 3"	435.00
50-1842	20 x 18 x 3"	462.50

Amara 250



- for moderate/high risk of skin breakdown; with foam insert
- moderate postural support

W x D x H

50-1810	16 x 16 x 3"	155.00
50-1811	18 x 18 x 3"	180.00
50-1812	20 x 18 x 3"	212.50

Amara solid seat inserts



- wood insert helps eliminate sagging and hammocking of wheelchair sling seats

W x D

50-1850	16 x 16"	35.00
50-1851	18 x 18"	37.50
50-1852	20 x 20"	40.00

high seat cushion



- use to elevate seating level of normal chair to reduce stress on knees and legs
- vinyl covered, wood base
- over shoulder carry strap for portability
- 16"W x 16"D x 4"H
- specify color: black (BLK), burgundy (BUR) or blue (B)

50-1880	high seat	60.00
---------	-----------	-------

foam cushions

- high-density foam for seating comfort
- includes removable, washable navy cover



W x D x H

50-1325	16 x 18 x 2"	20.00
50-1326	16 x 18 x 3"	22.50
50-1327	16 x 18 x 4"	25.00

gel cushions

- helps prevent and control decubitus problems
- gel for displacement of pressure
- removable navy cover




W x D x H

50-1360	16 x 18 x 2"	45.00
50-1361	16 x 20 x 2"	45.00

foam with gel bladder wheelchair cushions

- assists in prevention, treatment and management of pressure ulcers
- gel for pressure reduction, support and comfort
- fire retardant polyurethane foam shell
- removable/replaceable cover with water-resistant top and non-slip vinyl base
- 275 lb capacity



W x D x H

50-1362	16 x 18 x 2"	55.00
50-1364	16 x 16 x 3"	60.00
50-1365	16 x 18 x 3"	60.00
50-1366	16 x 20 x 3"	60.00