

bathroom products - toileting

elevated seat



- elevates seat by 4"
- wider than typical donut, slightly contoured
- hand grips for portability

43-2520	economy	25.00
43-2521	slip-in	35.00
43-2522	bolt-down	35.00

elevated seat with arms

- elevates seat by 4"
- for use with regular toilet seat
- width between arms is 19"; arm height is 4½"
- 300 lb capacity
- hardware included



43-2550	standard	57.50
43-2551	elongated	57.50

extra wide elevated seats with arms



- elevates seat by 4"
- arm rests and bolt-down bracket
- 22½" seat width; 18" width between arms
- for standard and elongated seats

- seat with aluminum legs for stability
- seat slides up legs to adjust height



A	43-2560	with arms	95.00
B	43-2565	with arms and legs	120.00

contoured elevated toilet seats

- provides for comfort and support
- available in three heights: 2", 4" and 6"
- standard seat includes a slip-in bracket for fast and safe installation
- elongated seat includes bolt-down bracket for a more permanent attachment



standard 6"

standard seat

43-2500	2" slip-in	42.50
43-2501	4" slip-in	45.00
43-2502	6" slip-in	47.50

elongated seat

43-2510	2" bolt-down	42.50
43-2511	4" bolt-down	45.00
43-2512	6" bolt-down	47.50

Lock-On™ seat with arms

- elevates seat by 5"
- openings in front and back for easy access
- width between arms is 19"; arm height is 4½"
- 300 lb capacity



43-2555	Lock-On™ seat w/arms	57.50
---------	----------------------	-------

Arthro Tall-ette seats



- elevates seat by 4"
- slopes on one side for those requiring extended leg positioning

right side sloped

43-2540R	standard slip-in	50.00
43-2541R	elongated bolt-down	50.00

left side sloped

43-2540L	standard slip-in	50.00
43-2541L	elongated bolt-down	50.00

splash guard



- extra-large guard directs urine flow towards toilet
- fits most standard and elongated toilets and elevated seats

45-1257	splash guard	12.50
---------	--------------	-------

hinged elevated toilet seats

- elevates seat by 4"
- two-part seat allows for normal raising and lowering of lid
- uses existing toilet seat and lid
- hinged seat allows for easy cleaning



43-2570	standard	67.50
43-2571	elongated	67.50

brackets

slip-in bracket

- slides under toilet seat hinge

43-2580	slip-in	12.50
---------	---------	-------



bolt-down bracket

- use toilet seat hinge bolts
- unscrew toilet seat hinge to slide in bracket

43-2581	bolt-down	15.00
---------	-----------	-------

