



lightweight cutlery

comfort grip cutlery

- lightweight, 3 oz
- soft vinyl coating with finger indentations for secure grip
- knife can also be used as a rocker knife and is suitable for one handed use



straight

61-0040	knife	10.00
61-0041	fork	10.00
61-0042	teaspoon	10.00
61-0043	soup spoon	10.00
61-0044	small teaspoon	10.00



bent specify left (L) or right (R)

61-0046	L/R fork	10.50
61-0047	L/R teaspoon	10.50
61-0048	L/R soup spoon	10.50



utensil holders

hand clip

- spring action hand clip fits snug against the hand
- attaches to objects with strips of heavy duty hook and loop
- use on phones, cups, bottles, hair brushes, etc.



holder with utensil pocket

- adjustable hook and loop webbing with utensil pocket



slide utensil into sleeve

61-0110	pediatric	7.00
61-0111	small/medium	7.00
61-0112	large	7.00

61-0100 hand clip 9.00



universal built-up handle

- ergonomic, soft, ribbed, built-up handle makes daily activities easier
- instantly creates larger gripping surface
- 1 3/4" (4.4 cm) diameter accommodates utensils with widths from 1/4"-1" (.64-2.54 cm)
- can also be used with everyday items such as pens, tooth brushes and knitting needles
- dishwasher safe



61-0008 universal handle, set of 4 20.00

soft handle cutlery

- large, foam grip handles
- lightweight, 1 oz

straight

61-0060	knife	13.50
61-0061	fork	13.50
61-0062	teaspoon	13.50
61-0063	soup spoon	13.50



bent

specify left (L) or right (R)

61-0066	L/R fork	14.50
61-0067	L/R teaspoon	14.50
61-0068	L/R soup spoon	14.50



lightweight cutlery

- lightweight, 1.7 oz
- 1" diameter handle is slightly tapered for easy grip
- knife can also be used as a rocker knife



straight

61-0050	knife	10.00
61-0051	fork	10.00
61-0052	teaspoon	10.00
61-0053	soup spoon	10.00

bent

specify left (L) or right (R)

61-0055	L/R knife	12.50
61-0056	L/R fork	11.00
61-0057	L/R teaspoon	11.00
61-0058	L/R soup spoon	11.00



utensils with finger loop

- has look and feel of traditional utensils
- designed for people with little to no grip strength



- user places finger or thumb through the loop and stabilizes the utensil with the palm

- vinyl-coated handles

- dishwasher safe

specify left (L) or right (R)

61-0004	L/R teaspoon	12.00
61-0005	L/R tablespoon	12.00
61-0006	L/R salad fork	12.00
61-0007	L/R dinner fork	12.00